

Sterilization Container Systems

Secure Today
Elevate Tomorrow







Different Lids Colors

Our sterilization containers come in eight different colors (silver, yellow, red, blue, green, black, bronze and purple) ensuring to easily distinguish between different sets or types of instruments at a glance.




Advantages of Multiple Colors

- **Enhanced Organization:** Different colors can be used to categorize and organize instruments based on their usage, department, or level of sterilization, reducing the risk of cross-contamination and improving efficiency..
- **Easy Identification:** Color-coding helps in quickly identifying the contents of each container, saving time and reducing errors in busy environments.
- **Improved Workflow:** With a clear visual system in place, staff can more easily manage inventory and ensure that the right instruments are always in the right place.

Customized Solutions & Tailored to Your Needs

Our sterilization containers are manufactured in 8 main sizes to accommodate various surgical needs: standard sizes (1/1 Full Size, 3/4 Quarter Size, 1/2 Half Size, Mini, Dental and 1/2 Dental Size) are ideal for common instrument sets, while the Wide Body size suits specialized instruments like Bookwalter retractors, and the Long Body size is designed for extended instruments such as robotic arms.

The sterilization containers are crafted in five primary heights to suit a variety of surgical instruments. Additionally, custom sizes or heights can be tailored to meet specific requirements.



Standard 8 Different Size Types
For Standard & Reusable Models

Type	Packaging according to EN 868-8		Steam sterilization procedure according to EN 285				
	Description	Symbolic Packaging	Gravity	Pre-Vacuum	Fractional Vacuum	Fractional Flow	Steam Injection
Sterilization Container with filters	Textiles in container with filter in lid and box		+	+	+	+	+
			O	O	+	+	+
	Textiles in container with filter in lid, non-perforated box		-	-	+	+	+
			-	-	+	+	O
	Textiles in container with filter in box, non-perforated lid		O	+	+	+	+
			-	O	+	+	+
	Instruments in wire basket, in sterilization filter in lid and box		+3	+	+	+	+
	Instruments in wire basket, in sterilization container filter in lid and non-perforated box		-	+2	+2	+2	+2
Sterilization Container with valves	Textiles in sterilization container with valves in lid and box		-	O	+	+	O
			-	-	+	+	O
	Textiles in sterilization container with valves in lid and no valves on the box		-	-	+	+	-
			-	-	+	+	-
	Instruments in wire basket in sterilization container with valve in lid and box		-	+2	+2	+2	+2
	Instruments in wire basket, in sterilization container with valve in lid and non-perforated box		-	+2	+2	+2	+2

*To prevent errors during service, operation, or loading of the autoclaves, consult above table.
*Use only sterilization containers that guarantee proper evaporation or discharge of condensate.
*Keep in mind that gravitation-based procedures, which lack a dedicated drying process, do not ensure sufficient after-drying effects.



STANDARD & REUSABLE MODELS

SECURE TODAY

ELEVATETOMORROW



WHY STERILIZATION CONTAINERS?

- Consistent Processing
- Reduced Risk of Contamination
- Physical Instrument Protection
- Optimal Set Organization
- Reusability / Sustainability
- Reduced Waste
- Streamlined Process
- Easy Identification and Access
- Adherence to Standards
- Minimized Handling

Today's Safety. Tomorrow's Efficiency



Elevated Sterilization Solutions: For Your Operational Excellence

Discover our 2mm special alloy aluminum container systems that provide superior protection for your surgical instruments. Available in 8 different color palettes, they facilitate organization and quick identification. This solution supports your operational excellence by combining reliability, an eco-friendly approach, and long-term cost advantages.

SAFETY LIDS



TWO DIFFERENT FILTER TYPES

Paper Filter or PTFE Filter



EASE WITH COLORS

EXCELLENCE IN WORKFLOW

 Different Lids Colors



TWO BOTTOM TYPES

 Non Perforated

 Perforated



Standard Model

1/1 Full Size

Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



Code Outside Dimension Inside Dimension L x W x H

K1-10-110 597 x 275 x 95 mm 546 x 258 x 90 mm


K1-10-135 597 x 275 x 120 mm 546 x 258 x 115 mm

K1-10-150 597 x 275 x 136 mm 546 x 258 x 131 mm

K1-10-200 597 x 275 x 188 mm 546 x 258 x 183 mm

K1-10-260 597 x 275 x 235 mm 546 x 258 x 230 mm

230 mm Lids



Colors	Code
Grey	K1-10-
Yellow	1 K1-
Red	10-2
Blue	K1-10-
Green	3 K1-
Black	10-4
Purple	K1-10-
Bronze	5 K1-
10-6	
Dimension	585 x 285 x 137 mm
7 K1-	
10-8	



Perforated



Code Outside Dimension Inside Dimension L x W x H

KP1-10-110 597 x 275 x 95 mm 544 x 257 x 75 mm

KP1-10-135 597 x 275 x 120 mm 544 x 257 x 100 mm

KP1-10-150 597 x 275 x 136 mm 544 x 257 x 116 mm

KP1-10-200 597 x 275 x 188 mm 544 x 257 x 168 mm

KP1-10-260 597 x 275 x 235 mm 544 x 257 x 215 mm

Safety Lids



Colors	Code
Grey	KS1-10-
Yellow	1 KS1-
Red	10-2
Blue	KS1-10-
Green	3 KS1-
Black	10-4
Purple	KS1-10-
Bronze	5 KS1-
10-6	
Dimension	586 x 289 x 136 mm
7 KS1-	
10-8	

Wire Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box


Code	Outside Dimension L x W x H
B1-5425-30	540 x 245 x 30 mm
B1-5425-50	540 x 245 x 50 mm
B1-5425-70	540 x 245 x 70 mm
B1-5425-100	540 x 245 x 100 mm




Lid

Code	Outside Dimension L x W x H
B1-5425	510 x 230 x 5 mm

Perforated Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box (with 10 mm feet)


Code	Outside Dimension L x W x H
B2-5425-30	540 x 251 x 40 mm
B2-5425-50	540 x 251 x 60 mm
B2-5425-70	540 x 251 x 80 mm
B2-5425-100	540 x 251 x 110 mm



Lid

Code	Outside Dimension L x W x H
B2-5425	540 x 256 x 5 mm

Hybrid Basket



Mesh: 4x4 mm Sides: ø 3 mm Wire: ø 1 mm



Box

Code	Outside Dimension L x W x H
B3-5425-50	540 x 255 x 50 mm
B3-5425-70	540 x 255 x 70 mm
B3-5425-100	540 x 255 x 100 mm



Lid

Code	Outside Dimension L x W x H
B3-5425	510 x 240 x 5 mm

Silicone Mat



Code	Outside Dimension L x W x H
A7-5425	500 x 235 mm


Standard Model

1/1 Full Size Baskets

Wire Basket



Perforated Basket



Hybrid Basket



Standard Model

3/4Quarter Size

 Different Lids
Colors



With Safety Lid



With Perforated Lid



Non
Perforated



Code Outside Dimension
 Inside Dimension
 L x W x H

K1-20-110 477 x 275 x 95 mm
 423 x 258 x 90 mm

K1-20-135 477 x 275 x 120 mm
 423 x 258 x 115 mm

K1-20-150 477 x 275 x 136 mm
 423 x 258 x 131 mm

K1-20-200 477 x 275 x 188 mm
 423 x 258 x 183 mm



Perforated




Code Outside Dimension
 Inside Dimension
 L x W x H

KP1-20-110 477 x 275 x 95 mm
 423 x 258 x 75 mm

KP1-20-135 477 x 275 x 120 mm
 423 x 258 x 100 mm

KP1-20-150 477 x 275 x 136 mm
 423 x 258 x 116 mm

KP1-20-200 477 x 275 x 188 mm
 423 x 258 x 168 mm



Lids



Colors Code

Grey K1-20-

Yellow 1 K1-

Red 20-2

Blue K1-20-

Green 3 K1-

Black 20-4

Purple K1-20-

Bronze 5 K1-

 20-6

 Dimension 465 x 289 x 137 mm

 7 K1-

 20-8



Safety Lids



Colors Code

Grey KS1-20-

Yellow 1 KS1-

Red 20-2

Blue KS1-20-

Green 3 KS1-

Black 20-4

Purple KS1-20-

Bronze 5 KS1-

 20-6

 Dimension 465 x 289 x 137 mm

 7 KS1-

 20-8

Wire Basket



Box

Code	Outside Dimension L x W x H
B1-4025-30	405 x 245 x 30 mm
B1-4025-50	405 x 245 x 50 mm
B1-4025-70	405 x 245 x 70 mm
B1-4025-100	405 x 245 x 100 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Lid

Code	Outside Dimension L x W x H
B1-4025	375 x 230 x 5 mm

Perforated Basket



Box (with 10 mm feet)

Code	Outside Dimension L x W x H
B2-4025-30	403 x 251 x 40 mm
B2-4025-50	403 x 251 x 60 mm
B2-4025-70	403 x 251 x 80 mm
B2-4025-100	403 x 251 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Lid

Code	Outside Dimension L x W x H
B2-4025	405 x 256 x 5 mm

Hybrid Basket



Box

Code	Outside Dimension L x W x H
B3-4025-50	405 x 255 x 50 mm
B3-4025-70	405 x 255 x 70 mm
B3-4025-100	405 x 255 x 100 mm


Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm



Lid

Code	Outside Dimension L x W x H
B3-4025	375 x 240 x 5 mm

Silicone Mat





Code	Outside Dimension L x W x H
A7-4025	370 x 235 mm

Standard Model

3/4 Quarter SizeBaskets


Wire Basket



Perforated Basket



Hybrid Basket



Standard Model

1/2 Half Size

Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



Code Outside Dimension Inside Dimension L x W x H

K1-30-110 307 x 275 x 95 mm 258 x 258 x 90 mm

K1-30-135 307 x 275 x 120 mm 258 x 258 x 115 mm

K1-30-150 307 x 275 x 136 mm 258 x 258 x 131 mm

K1-30-200 307 x 275 x 188 mm 258 x 258 x 183 mm

K1-30-260 307 x 275 x 235 mm 258 x 258 x 230 mm



Colors Code

Grey K1-30-

Yellow 1 K1-

Red 30-2

Blue K1-30-

Green 3 K1-

Black 30-4

Purple K1-30-

Bronze 5 K1-

30-6

Dimension 295 x 289 x 137 mm

7 K1-

30-8



Perforated



Code Outside Dimension Inside Dimension L x W x H

KP1-30-110 307 x 275 x 95 mm 258 x 258 x 75 mm

KP1-30-135 307 x 275 x 120 mm 258 x 258 x 100 mm

KP1-30-150 307 x 275 x 136 mm 258 x 258 x 116 mm

KP1-30-260 307 x 275 x 188 mm 258 x 258 x 168 mm

KP1-30-260 307 x 275 x 235 mm 258 x 258 x 215 mm



Colors Code

Grey KS1-30-

Yellow 1 KS1-

Red 30-2

Blue KS1-30-

Green 3 KS1-

Black 30-4

Purple KS1-30-

Bronze 5 KS1-


30-6

Dimension 295 x 289 x 137 mm

7 KS1-

30-8

Wire Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box

Code Outside Dimension L x W x H
B1-2525-30 246 x 246 x 30 mm
B1-2525-50 246 x 246 x 50 mm
B1-2525-70 246 x 246 x 70 mm
B1-2525-100 246 x 246 x 100 mm




Lid

Code B1-2525
Outside Dimension L x W x H 240 x 215 x 5 mm

Perforated Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box (with 10 mm feet)

Code Outside Dimension L x W x H
B2-2525-30 244 x 250 x 40 mm
B2-2525-50 244 x 250 x 60 mm
B2-2525-70 244 x 250 x 80 mm
B2-2525-100 244 x 250 x 110 mm




Lid

Code B2-2525
Outside Dimension L x W x H 247 x 256 x 5 mm

Hybrid Basket



Mesh: 4x4 mm Sides: ø 3 mm Wire: ø 1 mm



Box

Code Outside Dimension L x W x H
B3-2525-50 245 x 255 x 50 mm
B3-2525-70 245 x 255 x 70 mm
B3-2525-100 245 x 255 x 100 mm



Lid

Code B3-2525
Outside Dimension L x W x H 215 x 240 x 5 mm

Silicone Mat




Code A7-2525
Outside Dimension L x W x H 225 x 235 mm

Standard Model

1/2 Half Size Baskets


Wire Basket



Perforated Basket



Hybrid Basket



Standard Model

Long Size

Different Lids
Colors



With Perforated Lid



With Color Options



Non
Perforated



Code Outside Dimension
Inside Dimension
L x W x H

K1-40-150 770 x 275 x 136
mm 715 x 258 x
131 mm

K1-40-200 770 x 275 x 188
mm 715 x 258 x
183 mm



Perforated



Code Outside Dimension
Inside Dimension
L x W x H

KP1-40-150 770 x 275 x 136
mm 715 x 258 x
116 mm

KP1-40-200 770 x 275 x 188
mm 715 x 258 x
168 mm

Lids

Colors Code

Grey K1-40-1

Yellow K1-40-2

Red K1-40-3

Blue K1-40-4

Green K1-40-5

Black K1-40-6

Purple K1-40-7

Bronze K1-40-8

Dimension
758 x 280 x 37 mm



Wire Basket



Box

Code Outside Dimension
L x W x H

B1-6625-50 665 x 245 x 50 mm

B1-6625-70 665 x 245 x 70 mm

B1-6625-100 665 x 245 x 100 mm


Mesh: 4x4 mm
Wire: ø 1 mm

Lid

Code
B1-6625

Outside Dimension
L x W x H
635 x 230 x 5 mm

Perforated Basket



Box (with 10 mm feet)

Code Outside Dimension
L x W x H

B2-6625-30 650 x 251 x 40 mm

B2-6625-50 650 x 251 x 60 mm

B2-6625-70 650 x 251 x 80 mm

B2-6625-100 650 x 251 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm

Lid

Code
B2-6625

Outside Dimension
L x W x H
652 x 256 x 5 mm

Hybrid Basket



Box

Code Outside Dimension
L x W x H

B3-6625-50 665 x 255 x 50 mm

B3-6625-70 665 x 255 x 70 mm

B3-6625-100 665 x 255 x 100 mm

Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm

Lid


Code
B3-6625

Outside Dimension
L x W x H
635 x 240 x 5 mm


Standard Model

Long Size Baskets

Wire Basket



Perforated Basket



Hybrid Basket



Standard Model

Wide Size

 Different Lids
Colors



With Perforated Lid



With Color Options



Non
Perforated



Code Outside Dimension
 Inside Dimension
 L x W x H

K1-50-110 597 x 365 x 105
 mm 546 x 348 x
 100 mm
K1-50-170 597 x 365 x 165
 mm 546 x 348 x
 160 mm




Perforated



Code Outside Dimension
 Inside Dimension
 L x W x H

KP1-50-110 597 x 365 x 105
 mm 546 x 348 x 85
 mm
KP1-50-170 597 x 365 x 165
 mm 546 x 348 x
 145 mm

 Lids

Colors Code

Grey K1-50-1

Yellow K1-50-2

Red K1-50-3

Blue K1-50-4

Green K1-50-5


Black K1-50-6

Purple K1-50-7

Bronze K1-50-8

Dimension
585 x 370 x 37 mm

Wire Basket



 Box


Code	Outside Dimension L x W x H
B1-5333-50	533 x 337 x 50 mm
B1-5333-70	533 x 337 x 70 mm
B1-5333-100	533 x 337 x 100 mm

Mesh: 4x4 mm
Wire: ø 1 mm

 Lid

Code
B1-5333
Outside Dimension
L x W x H
503 x 322 x 5 mm

Perforated Basket



 Box (with 10 mm feet)


Code	Outside Dimension L x W x H
B2-5333-30	535 x 340 x 40 mm
B2-5333-50	535 x 340 x 60 mm
B2-5333-70	535 x 340 x 80 mm
B2-5333-100	535 x 340 x 110 mm


Mesh: 4x4 mm
Wire: ø 1 mm

 Lid

Code
B2-5333
Outside Dimension
L x W x H
537 x 343 x 5 mm

Hybrid Basket



 Box

Code	Outside Dimension L x W x H
B3-5333-50	533 x 337 x 50 mm
B3-5333-70	533 x 337 x 70 mm
B3-5333-100	533 x 337 x 100 mm

Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm


 Lid

Code
B3-5333
Outside Dimension
L x W x H
503 x 322 x 5 mm

Standard Model

Wide Size Baskets

Wire Basket



Perforated Basket



Hybrid Basket



Standard Model

Flat Size

 Different Lids
Colors



With Perforated Lid



With Color Options



Non
Perforated

Code	Outside Dimension Inside Dimension L x W x H
K1-90-45	297 x 275 x 55 mm 258 x 258 x 50 mm
K1-90-75	297 x 275 x 75 mm 258 x 258 x 70 mm




Lids

Colors	Code
Grey	K1-90-1
Yellow	K1-90-2
Red	K1-90-3
Blue	K1-90-4
Green	K1-90-5
Black	K1-90-6
Purple	K1-90-7
Bronze	K1-90-8

Dimension
295 x 280 x 37 mm

Wire Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Box

Code	Outside Dimension L x W x H
B1-2525-30	246 x 246 x 30 mm
B1-2525-50	246 x 246 x 50 mm



Lid

Code	Outside Dimension L x W x H
B1-2525	240 x 215 x 5 mm

Perforated Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Box (with 10 mm feet)

Code	Outside Dimension L x W x H
B2-2525-30	244 x 250 x 40 mm
B2-2525-50	244 x 250 x 60 mm



Lid

Code	Outside Dimension L x W x H
B2-2525	247 x 256 x 5 mm

Silicone Mat




Code	Outside Dimension L x W x H
A7-2525	225 x 235 mm

Standard Model


Flat Size Baskets



Wire Basket



Perforated Basket



Standard Model

MiniSize

 Different Lids
Colors



With Perforated Lid




With Color Options



Non
Perforated

Code	Outside Dimension Inside Dimension L x W x H
K1-70-40	305 x 135 x 45 mm 272 x 125 x 33 mm
K1-70-70	305 x 135 x 60 mm 272 x 125 x 48 mm
K1-70-100	305 x 135 x 90 mm 272 x 125 x 78 mm
K1-70-130	305 x 135 x 120 mm 272 x 125 x 108 mm



Lids


Colors	Code
Grey	K1-70-1
Yellow	K1-70-2
Red	K1-70-3
Blue	K1-70-4
Green	K1-70-5
Black	K1-70-6
Purple	K1-70-7
Bronze	K1-70-8

Dimension
303 x 136 x 29 mm

Wire Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Box

Code	Outside Dimension L x W x H
B1-2712-40	268 x 117 x 40 mm
B1-2712-100	268 x 117 x 100 mm



Lid

Code	Outside Dimension L x W x H
B1-2712	250 x 120 x 5 mm

Perforated Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Box (with 10 mm feet)

Code	Outside Dimension L x W x H
B2-2712-40	260 x 120 x 50 mm
B2-2712-100	260 x 120 x 110 mm



Lid

Code	Outside Dimension L x W x H
B2-2712	262 x 122 x 5 mm

Silicone Mat



Code	Outside Dimension L x W x H
A7-2712	270 x 125 mm


Standard Model

MiniSizeBaskets

Wire Basket



Perforated Basket



Standard Model

Dental Size

 Different Lids
Colors



With Perforated Lid



With Color Options




Code	Outside Dimension Inside Dimension L x W x H
K1-60-40	310 x 189 x 45 mm 277 x 180 x 33 mm
K1-60-70	310 x 189 x 60 mm 277 x 180 x 48 mm
K1-60-100	310 x 189 x 90 mm 277 x 180 x 78 mm
K1-60-130	310 x 189 x 120 mm 277 x 180 x 108 mm



Colors	Code
Grey	K1-60-1
Yellow	K1-60-2
Red	K1-60-3
Blue	K1-60-4
Green	K1-60-5
Black	K1-60-6
Purple	K1-60-7
Bronze	K1-60-8

Dimension
308 x 190 x 29 mm

Wire Basket




Code	Outside Dimension L x W x H
B1-2717-40	270 x 170 x 40 mm
B1-2717-100	270 x 170 x 100 mm

Mesh: ø 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B1-2717	252 x 172 x 5 mm

Perforated Basket



Code	Outside Dimension L x W x H
B2-2717-40	273 x 172 x 50 mm
B2-2717-100	273 x 172 x 110 mm

Mesh: ø 4x4 mm
Wire: ø 1 mm




Code	Outside Dimension L x W x H
B2-2717	275 x 174 x 5 mm

Standard Model

Dental Size Baskets

Wire Basket



Perforated Basket



Standard Model

1/2 Dental Size

 Different Lids
Colors



With Perforated Lid



With Color Options



Non
Perforated

Code	Outside Dimension Inside Dimension L x W x H
K1-80-40	207 x 144 x 45 mm 174 x 135 x 33 mm
K1-80-70	207 x 144 x 60 mm 174 x 135 x 48 mm
K1-80-100	207 x 144 x 90 mm 174 x 135 x 78 mm
K1-80-130	207 x 144 x 120 mm 174 x 135 x 108 mm



Lids

Colors	Code
Grey	K1-80-1
Yellow	K1-80-2
Red	K1-80-3
Blue	K1-80-4
Green	K1-80-5
Black	K1-80-6
Purple	K1-80-7
Bronze	K1-80-8

Dimension
205 x 145 x 29 mm

Wire Basket



Mesh: $\varnothing$ 4x4 mm
Wire: $\varnothing$ 1 mm



Box

Code	Outside Dimension L x W x H
B1-1714-40	165 x 130 x 35 mm
B1-1714-100	165 x 130 x 100 mm



Lid

Code	Outside Dimension L x W x H
B1-1714	147 x 132 x 5 mm

Perforated Basket



Mesh: $\varnothing$ 4x4 mm
Wire: $\varnothing$ 1 mm



Box (with 10 mm feet)

Code	Outside Dimension L x W x H
B2-1714-40	165 x 130 x 45 mm
B2-1714-100	165 x 130 x 100 mm



Lid

Code	Outside Dimension L x W x H
B2-1714	167 x 130 x 5 mm

Standard Model

1/2 Dental Size Baskets

Wire Basket




Perforated Basket



Clean. Safe. Sustainable.

2 Different Filtering System
Mechanical Valve System & Barrier System




Different Lids Colors

Danhealth.com

BARRIER MODEL

The Barrier Model is an innovative sterilization container designed to ensure maximum safety in CSSD workflows. Its labyrinth structure prevents microorganisms from entering, maintaining sterility for extended periods. This design offers reliable, practical, and scientifically proven protection.


BARRIER FILTERS



EASY-TO-USE LOCK / UNLOCK



STEAM FILTER CROSS-SECTIONAL VIEW



EASE WITH COLORS
EXCELLENCE IN WORKFLOW

Different Lids Colors



Scientific Foundation

The principle behind the Barrier mechanism traces back to Louis Pasteur's groundbreaking "Swan Neck Flask" experiment in 1862, which challenged the idea of spontaneous generation. Pasteur designed a flask with an S-shaped neck that allowed air to pass freely but trapped dust and microorganisms in its curves. When the broth inside was boiled and left untouched, no microbial growth occurred, clearly proving that life could not emerge without contamination from existing microorganisms. However, once the flask was tilted and the trapped particles reached the broth, microbial growth appeared almost immediately.

The Barrier Model adopts the same principle through its labyrinth-like system. By forcing air to travel through a complex path, it prevents microorganisms and viruses from penetrating the container. This ensures that sterility is maintained, while still allowing necessary airflow. The result is a reliable and scientifically validated protective system, directly inspired by one of the most important experiments in microbiology.

REDUCE. REUSE. RENEW.

Seize Cost Control and Sustainability: Ditch the waste expense and embrace the smart, permanent solution.



VALVE MODEL

Valve Model sterilization containers are designed with an advanced mechanical valve system that removes the need for disposable filters while delivering stable and reliable performance throughout extended use.

Each valve is precisely engineered with integrated gaskets that ensure smooth operation in response to pressure variations during the sterilization process. The secure closure, supported by high-grade silicone gaskets, provides excellent sealing integrity and long-term durability.


Combining innovative design with practicality, Dans Valve Model containers offer a sustainable and cost-efficient solution for safe and effective sterilization in healthcare environments worldwide.

EASE WITH COLORS EXCELLENCE IN WORKFLOW

 Different Lids Colors



RUNNING PRINCIPLE



At the beginning of the sterilization cycle, the vacuum phase causes the valves to fully open, allowing the evacuation of dry air from the container.

During the sterilization cycle, the pressure-activated valves automatically open under positive pressure, enabling controlled steam penetration for effective sterilization.

Once the process is complete, and no external forces such as vacuum or pressure are present, the entire system remains sealed, preserving the sterility of its contents.

Reusable Model

1/1 Full Size

 Different Lids
Colors



With Reusable Mechanism




With Color Options




Non
Perforated

Code	Outside Dimension Inside Dimension L x W x H
K1-10-110	597 x 275 x 95 mm 546 x 258 x 90 mm
K1-10-135	597 x 275 x 120 mm 546 x 258 x 115 mm
K1-10-150	597 x 275 x 136 mm 546 x 258 x 131 mm
K1-10-200	597 x 275 x 188 mm 546 x 258 x 183 mm
K1-10-260	597 x 275 x 235 mm 546 x 258 x 230 mm



Valve Lids

Colors	Code
Grey	K2-10-
Yellow	1 K2-
Red	10-2
Blue	K2-10-
Green	3 K2-
Black	10-4
Purple	K2-10-
Bronze	5 K2-
Dimension	10-6
585 x 280 x 210 mm	7 K2-
	10-8



Barrier Lids

Colors	Code
Grey	K3-10-
Yellow	1 K3-
Red	10-2
Blue	K3-10-
Green	3 K3-
Black	10-4
Purple	K3-10-
Bronze	5 K3-
Dimension	10-6
585 x 280 x 210 mm	7 K3-
	10-8

Wire Basket



Box

Code	Outside Dimension L x W x H
B1-5425-30	540 x 245 x 30 mm
B1-5425-50	540 x 245 x 50 mm
B1-5425-70	540 x 245 x 70 mm
B1-5425-100	540 x 245 x 100 mm


Mesh: 4x4 mm
Wire: ø 1 mm



Lid

Code	Outside Dimension L x W x H
B1-5425	510 x 230 x 5 mm

Perforated Basket



Box (with 10 mm feet)

Code	Outside Dimension L x W x H
B2-5425-30	540 x 251 x 40 mm
B2-5425-50	540 x 251 x 60 mm
B2-5425-70	540 x 251 x 80 mm
B2-5425-100	540 x 251 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Lid

Code	Outside Dimension L x W x H
B2-5425	540 x 256 x 5 mm

Hybrid Basket



Box

Code	Outside Dimension L x W x H
B3-5425-50	540 x 255 x 50 mm
B3-5425-70	540 x 255 x 70 mm
B3-5425-100	540 x 255 x 100 mm


Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm



Lid

Code	Outside Dimension L x W x H
B3-5425	510 x 240 x 5 mm

Silicone Mat




Code	Outside Dimension L x W x H
A7-5425	500 x 235 mm

Reusable Model

1/1 Full Size Baskets


Wire Basket



Perforated Basket



Hybrid Basket



Reusable Model

3/4 Quarter Size

Different Lids
Colors



With Reusable Mechanism



With Color Options



Code	Outside Dimension Inside Dimension L x W x H
K1-20-110	477 x 275 x 95 mm 423 x 258 x 90 mm
K1-20-135	477 x 275 x 120 mm 423 x 258 x 115 mm
K1-20-150	477 x 275 x 136 mm 423 x 258 x 131 mm
K1-20-200	477 x 275 x 188 mm 423 x 258 x 183 mm



Colors	Code
Grey	K2-20-
Yellow	1 K2-
Red	20-2
Blue	K2-20-
Green	3 K2-
Black	20-4
Purple	K2-20-
Bronze	5 K2-
Dimension	20-6
465 x 280 x 20 mm	7 K2-
	20-8



Colors	Code
Grey	K3-20-
Yellow	1 K3-
Red	20-2
Blue	K3-20-
Green	3 K3-
Black	20-4
Purple	K3-20-
Bronze	5 K3-
Dimension	20-6
465 x 280 x 20 mm	7 K3-
	20-8

Wire Basket



Code	Outside Dimension L x W x H
B1-4025-30	405 x 245 x 30 mm
B1-4025-50	405 x 245 x 50 mm
B1-4025-70	405 x 245 x 70 mm
B1-4025-100	405 x 245 x 100 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B1-4025	375 x 230 x 5 mm

Perforated Basket




Code	Outside Dimension L x W x H
B2-4025-30	403 x 251 x 40 mm
B2-4025-50	403 x 251 x 60 mm
B2-4025-70	403 x 251 x 80 mm
B2-4025-100	403 x 251 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B2-4025	405 x 256 x 5 mm

Hybrid Basket




Code	Outside Dimension L x W x H
B3-4025-50	405 x 255 x 50 mm
B3-4025-70	405 x 255 x 70 mm
B3-4025-100	405 x 255 x 100 mm

Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B3-4025	375 x 240 x 5 mm

Silicone Mat




Code	Outside Dimension L x W x H
A7-4025	370 x 235 mm


Reusable Model

3/4 Quarter Size Baskets


Wire Basket



Perforated Basket



Hybrid Basket



Reusable Model

1/2 Half Size

 Different Lids
Colors



With Reusable Mechanism




With Color Options



Non
Perforated

Code	Outside Dimension Inside Dimension L x W x H
K1-30-110	307 x 275 x 95 mm 258 x 258 x 90 mm
K1-30-135	307 x 275 x 120 mm 258 x 258 x 115 mm
K1-30-150	307 x 275 x 136 mm 258 x 258 x 131 mm
K1-30-200	307 x 275 x 188 mm 258 x 258 x 183 mm
K1-30-260	307 x 275 x 235 mm 258 x 258 x 230 mm



Valve Lids

Colors	Code
Grey	K2-30-
Yellow	1 K2-
Red	30-2
Blue	K2-30-
Green	3 K2-
Black	30-4
Purple	K2-30-
Bronze	5 K2-
Dimension	30-6
295 x 280 x 230 mm	



Barrier Lids

Colors	Code
Grey	K3-30-
Yellow	1 K3-
Red	30-2
Blue	K3-30-
Green	3 K3-
Black	30-4
Purple	K3-30-
Bronze	5 K3-
Dimension	30-6
295 x 280 x 230 mm	
	7 K3-
	30-8

Wire Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Box


Code	Outside Dimension L x W x H
B1-2525-30	246 x 246 x 30 mm
B1-2525-50	246 x 246 x 50 mm
B1-2525-70	246 x 246 x 70 mm
B1-2525-100	246 x 246 x 100 mm



Lid

Code	Outside Dimension L x W x H
B1-2525	240 x 215 x 5 mm

Perforated Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Box (with 10 mm feet)


Code	Outside Dimension L x W x H
B2-2525-30	244 x 250 x 40 mm
B2-2525-50	244 x 250 x 60 mm
B2-2525-70	244 x 250 x 80 mm
B2-2525-100	244 x 250 x 110 mm



Lid

Code	Outside Dimension L x W x H
B2-2525	247 x 256 x 5 mm

Hybrid Basket



Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm



Box

Code	Outside Dimension L x W x H
B3-2525-50	245 x 255 x 50 mm
B3-2525-70	245 x 255 x 70 mm
B3-2525-100	245 x 255 x 100 mm



Lid

Code	Outside Dimension L x W x H
B3-2525	215 x 240 x 5 mm

Silicone Mat




Code	Outside Dimension L x W x H
A7-2525	225 x 235 mm

Reusable Model

1/2 Half Size Baskets


Wire Basket



Perforated Basket



Hybrid Basket



Reusable Model

Long Size

 Different Lids
Colors



With Reusable Mechanism



With Color Options



Code	Outside Dimension Inside Dimension L x W x H
K1-40-150	770 x 275 x 136 mm 715 x 258 x 131 mm
K1-40-200	770 x 275 x 188 mm 715 x 258 x 183 mm



Colors	Code
Grey	K2-40-
Yellow	1 K2-
Red	40-2
Blue	K2-40-
Green	3 K2-
Black	40-4
Purple	K2-40-
Bronze	5 K2-
Dimension	40-6 758 x 280 x 37 mm
	7 K2-
	40-8



Colors	Code
Grey	K3-40-
Yellow	1 K3-
Red	40-2
Blue	K3-40-
Green	3 K3-
Black	40-4
Purple	K3-40-
Bronze	5
Dimension	K3-40- 758 x 280 x 37 mm
	K3-40-
	7
	K3-40-
	8

Wire Basket



Code	Outside Dimension L x W x H
B1-6625-50	665 x 245 x 50 mm
B1-6625-70	665 x 245 x 70 mm
B1-6625-100	665 x 245 x 100 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B1-6625	635 x 230 x 5 mm

Perforated Basket



Code	Outside Dimension L x W x H
B2-6625-30	650 x 251 x 40 mm
B2-6625-50	650 x 251 x 60 mm
B2-6625-70	650 x 251 x 80 mm
B2-6625-100	650 x 251 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B2-6625	652 x 256 x 5 mm

Hybrid Basket



Code	Outside Dimension L x W x H
B3-6625-50	665 x 255 x 50 mm
B3-6625-70	665 x 255 x 70 mm
B3-6625-100	665 x 255 x 100 mm

Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm




Code	Outside Dimension L x W x H
B3-6625	635 x 240 x 5 mm


Reusable Model

Long Size Baskets

Wire Basket



Perforated Basket



Hybrid Basket



Reusable Model

Wide Size


Different Lids Colors



With Reusable Mechanism




With Color Options



Code	Outside Dimension Inside Dimension L x W x H
K1-50-110	597 x 365 x 105 mm 546 x 348 x 100 mm
K1-50-170	597 x 365 x 165 mm 546 x 348 x 160 mm



Colors	Code
Grey	K2-50-
Yellow	1 K2-
Red	50-2
Blue	K2-50-
Green	3 K2-
Black	50-4
Purple	K2-50-
Bronze	5 K2-
Dimension	50-6
585 x 370	K2-50-
7	K2-
50-8	



Colors	Code
Grey	K3-50-
Yellow	1 K3-
Red	50-2
Blue	K3-50-
Green	3 K3-
Black	50-4
Purple	K3-50-
Bronze	5 K3-
Dimension	50-6
585 x 370	K3-50-
7	K3-
50-8	

Wire Basket



Code	Outside Dimension L x W x H
B1-5333-50	533 x 337 x 50 mm
B1-5333-70	533 x 337 x 70 mm
B1-5333-100	533 x 337 x 100 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B1-5333	503 x 322 x 5 mm

Perforated Basket



Code	Outside Dimension L x W x H
B2-5333-30	535 x 340 x 40 mm
B2-5333-50	535 x 340 x 60 mm
B2-5333-70	535 x 340 x 80 mm
B2-5333-100	535 x 340 x 110 mm

Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B2-5333	537 x 343 x 5 mm

Hybrid Basket



Code	Outside Dimension L x W x H
B3-5333-50	533 x 337 x 50 mm
B3-5333-70	533 x 337 x 70 mm
B3-5333-100	533 x 337 x 100 mm

Mesh: 4x4 mm
Sides: ø 3 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B3-5333	503 x 322 x 5 mm

Reusable Model

Wide Size Baskets

Wire Basket



Perforated Basket



Hybrid Basket



Reusable Model

Flat Size

 Different Lids
Colors



With Reusable Mechanism



With Color Options



Code	Outside Dimension Inside Dimension L x W x H
K1-90-75	297 x 275 x 75 mm 258 x 258 x 70 mm



Colors	Code
Grey	K3-90-1
Yellow	K3-90-2
Red	K3-90-3
Blue	K3-90-4
Green	K3-90-5
Black	K3-90-6
Purple	K3-90-7
Bronze	K3-90-8

Dimension
295 x 280 x 37 mm

Wire Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B1-2525-30	246 x 246 x 30 mm
B1-2525-50	246 x 246 x 50 mm




Code	Outside Dimension L x W x H
B1-2525	240 x 215 x 5 mm

Perforated Basket



Mesh: 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B2-2525-30	244 x 250 x 40 mm
B2-2525-50	244 x 250 x 60 mm



Code	Outside Dimension L x W x H
B2-2525	247 x 256 x 5 mm


Silicone Mat




Code	Outside Dimension L x W x H
A7-2525	225 x 235 mm

Reusable Model

Flat Size Baskets



Wire Basket



Perforated Basket


Reusable Model

Mini Size

 Different Lids
Colors




With Reusable Mechanism



With Color Options



	Non Perforated
Code	Outside Dimension Inside Dimension L x W x H
K1-70-70	305 x 135 x 60 mm 272 x 125 x 48 mm
K1-70-100	305 x 135 x 90 mm 272 x 125 x 78 mm
K1-70-130	305 x 135 x 120 mm 272 x 125 x 108 mm

 Lids

Colors	Code
Grey	K3-70-1
Yellow	K3-70-2
Red	K3-70-3
Blue	K3-70-4
Green	K3-70-5
Black	K3-70-6
Purple	K3-70-7
Bronze	K3-70-8
Dimension	303 x 136 x 29 mm

Wire Basket



Mesh: 4x4 mm
Wire: ø 1 mm

 Box

 Lid

Code	Outside Dimension L x W x H
B1-2712-40	268 x 117 x 40 mm
B1-2712-100	268 x 117 x 100 mm

Code	Outside Dimension L x W x H
B1-2712	250 x 120 x 5 mm

Perforated Basket



Mesh: 4x4 mm
Wire: ø 1 mm

 Box (with 10 mm feet)

 Lid

Code	Outside Dimension L x W x H
B2-2712-40	260 x 120 x 50 mm
B2-2712-100	260 x 120 x 110 mm

Code	Outside Dimension L x W x H
B2-2712	262 x 122 x 5 mm

Silicone Mat



Code	Outside Dimension L x W x H
A7-2712	270 x 125 mm


Reusable Model

Mini Size Baskets

Wire Basket



Perforated Basket



Reusable Model

Dental Size


 Different Lids
Colors



With Reusable Mechanism



With Color Options



Code Outside Dimension
 Inside Dimension
 L x W x H


K1-60-70	310 x 189 x 60 mm 277 x 180 x 48 mm
K1-60-100	310 x 189 x 90 mm 277 x 180 x 78 mm
K1-60-130	310 x 189 x 120 mm 277 x 180 x 108 mm



Colors	Code
Grey	K3-60-1
Yellow	K3-60-2
Red	K3-60-3
Blue	K3-60-4
Green	K3-60-5
Black	K3-60-6
Purple	K3-60-7
Bronze	K3-60-8

Dimension
308 x 190 x 29 mm

Wire Basket




Code	Outside Dimension L x W x H
B1-2717-40	270 x 170 x 40 mm
B1-2717-100	270 x 170 x 100 mm

Mesh: ø 4x4 mm
Wire: ø 1 mm



Code	Outside Dimension L x W x H
B1-2717	252 x 172 x 5 mm

Perforated Basket



Code	Outside Dimension L x W x H
B2-2717-40	273 x 172 x 50 mm
B2-2717-100	273 x 172 x 110 mm

Mesh: ø 4x4 mm
Wire: ø 1 mm




Code	Outside Dimension L x W x H
B2-2717	275 x 174 x 5 mm

Reusable Model

Dental Size Baskets

Wire Basket



Perforated Basket



Reusable Model

1/2 Dental Size

 Different Lids
Colors



With Reusable Mechanism




With Color Options



Non
Perforated

Code	Outside Dimension Inside Dimension L x W x H
K1-80-70	207 x 144 x 60 mm 174 x 135 x 48 mm
K1-80-100	207 x 144 x 90 mm 174 x 135 x 78 mm
K1-80-130	207 x 144 x 120 mm 174 x 135 x 108 mm

 Lids

Colors	Code
Grey	K3-80-1
Yellow	K3-80-2
Red	K3-80-3
Blue	K3-80-4
Green	K3-80-5
Black	K3-80-6
Purple	K3-80-7
Bronze	K3-80-8

Dimension
205 x 145 x 29 mm

Wire Basket



Mesh: ø 4x4 mm
Wire: ø 1 mm



Box



Lid

Code	Outside Dimension L x W x H
B1-1714-40	165 x 130 x 35 mm
B1-1714-100	165 x 130 x 100 mm

Code	Outside Dimension L x W x H
B1-1714	147 x 132 x 5 mm

Perforated Basket



Mesh: ø 4x4 mm
Wire: ø 1 mm



Box (with 10 mm feet)



Lid

Code	Outside Dimension L x W x H
B2-1714-40	165 x 130 x 45 mm
B2-1714-100	165 x 130 x 110 mm

Code	Outside Dimension L x W x H
B2-1714	167 x 130 x 5 mm

Reusable Model

1/2 Dental Size Baskets

Wire Basket




Perforated Basket



Organize.
Sterilize.

Excel.

Holding the
Standard




Three Different Types of Basket
Models by Construction/Material



Perforated
Sheet

Ø 1 mm


□ 4 mm



Wire Mesh

Ø 1 mm

□ 4 mm



Hybrith

Ø 1 mm

□ 4 mm

Side Mesh

Ø 3 mm




Half Size




Box Code	Outside Dimension L x W x H
B1-2525-30	246 x 246 x 30 mm
B1-2525-50	246 x 246 x 50 mm
B1-2525-70	246 x 246 x 70 mm
B1-2525-100	246 x 246 x 100 mm
Lid Code	Outside Dimension L x W x H
B1-2525	240 x 215 x 5 mm

Long Size




Box Code	Outside Dimension L x W x H
B1-6625-50	665 x 245 x 50 mm
B1-6625-70	665 x 245 x 70 mm
B1-6625-100	665 x 245 x 100 mm
Lid Code	Outside Dimension L x W x H
B1-6625	635 x 230 x 5 mm

Wide Size




Box Code	Outside Dimension L x W x H
B1-5333-50	533 x 337 x 50 mm
B1-5333-70	533 x 337 x 70 mm
B1-5333-100	533 x 337 x 100 mm
Lid Code	Outside Dimension L x W x H
B1-5333	503 x 322 x 5 mm

Full Size




Box Code	Outside Dimension L x W x H
B1-5425-30	540 x 245 x 30 mm
B1-5425-50	540 x 245 x 50 mm
B1-5425-70	540 x 245 x 70 mm
B1-5425-100	540 x 245 x 100 mm
Lid Code	Outside Dimension L x W x H
B1-5425	510 x 230 x 5 mm

Din Size




Box Code	Outside Dimension L x W x H
B1-4825-30	480 x 245 x 30 mm
B1-4825-50	480 x 245 x 50 mm
B1-4825-70	480 x 245 x 70 mm
B1-4825-100	480 x 245 x 100 mm
Lid Code	Outside Dimension L x W x H
B1-4825	450 x 230 x 5 mm

Quarter Size



Box Code	Outside Dimension L x W x H
B1-4025-30	405 x 245 x 30 mm
B1-4025-50	405 x 245 x 50 mm
B1-4025-70	405 x 245 x 70 mm
B1-4025-100	405 x 245 x 100 mm
Lid Code	Outside Dimension L x W x H
B1-4025	375 x 230 x 5 mm

Dental Size



Box Code	Outside Dimension L x W x H
B1-2717-40	270 x 170 x 40 mm
B1-2717-100	270 x 170 x 100 mm
-	-
-	-
Lid Code	Outside Dimension L x W x H
B1-2717	252 x 172 x 5 mm

Mini Size




Box Code	Outside Dimension L x W x H
B1-2712-40	268 x 117 x 40 mm
B1-2712-100	268 x 117 x 100 mm
-	-
-	-
Lid Code	Outside Dimension L x W x H
B1-2712	250 x 120 x 5 mm

Dental 1/2



Box Code	Outside Dimension L x W x H
B1-1714-40	165 x 130 x 35 mm
B1-1714-100	165 x 130 x 100 mm
-	-
-	-
Lid Code	Outside Dimension L x W x H
B1-1714	147 x 132 x 5 mm



Half Size



Box Code	Outside Dimension L x W x H
B2-2525-30	244 x 250 x 40 mm
B2-2525-50	244 x 250 x 60 mm
B2-2525-70	244 x 250 x 80 mm
B2-2525-100	244 x 250 x 110 mm
Lid Code	Outside Dimension L x W x H
B2-2525	247 x 256 x 5 mm

Long Size



Box Code	Outside Dimension L x W x H
B2-6625-30	650 x 251 x 40 mm
B2-6625-50	650 x 251 x 60 mm
B2-6625-70	650 x 251 x 80 mm
B2-6625-100	650 x 251 x 110 mm
Lid Code	Outside Dimension L x W x H
B2-6625	652 x 256 x 5 mm

Wide Size



Box Code	Outside Dimension L x W x H
B2-5333-30	535 x 340 x 40 mm
B2-5333-50	535 x 340 x 60 mm
B2-5333-70	535 x 340 x 80 mm
B2-5333-100	535 x 340 x 110 mm
Lid Code	Outside Dimension L x W x H
B2-5333	537 x 343 x 5 mm

Full Size



Box Code	Outside Dimension L x W x H
B2-5425-30	540 x 251 x 40 mm
B2-5425-50	540 x 251 x 60 mm
B2-5425-70	540 x 251 x 80 mm
B2-5425-100	540 x 251 x 110 mm
Lid Code	Outside Dimension L x W x H
B2-5425	540 x 256 x 5 mm

Din Size



Box Code	Outside Dimension L x W x H
B2-4825-30	483 x 251 x 40 mm
B2-4825-50	483 x 251 x 60 mm
B2-4825-70	483 x 251 x 80 mm
B2-4825-100	483 x 251 x 110 mm
Lid Code	Outside Dimension L x W x H
B2-4825	485 x 256 x 5 mm

Quarter Size



Box Code	Outside Dimension L x W x H
B2-4025-30	403 x 251 x 40 mm
B2-4025-50	403 x 251 x 60 mm
B2-4025-70	403 x 251 x 80 mm
B2-4025-100	403 x 251 x 110 mm
Lid Code	Outside Dimension L x W x H
B2-4025	405 x 256 x 5 mm

Dental Size



Box Code	Outside Dimension L x W x H
B2-2717-40	273 x 172 x 50 mm
B2-2717-100	273 x 172 x 110 mm
-	-
-	-
Lid Code	Outside Dimension L x W x H
B2-2717	275 x 174 x 5 mm

Mini Size




Box Code	Outside Dimension L x W x H
B2-2712-40	260 x 120 x 50 mm
B2-2712-100	260 x 120 x 110 mm
-	-
-	-
Lid Code	Outside Dimension L x W x H
B2-2712	262 x 122 x 5 mm

Dental 1/2



Box Code	Outside Dimension L x W x H
B2-1714-40	165 x 130 x 45 mm
B2-1714-100	165 x 130 x 110 mm
-	-
-	-
Lid Code	Outside Dimension L x W x H
B2-1714	167 x 130 x 5 mm



Half Size

Box Code	Outside Dimension L x W x H
B3-2525-50	245 x 255 x 50 mm
B3-2525-70	245 x 255 x 70 mm
B3-2525-100	245 x 255 x 100 mm

Lid Code	Outside Dimension L x W x H
B3-2525	215 x 240 x 5 mm

Long Size

Box Code	Outside Dimension L x W x H
B3-6625-50	665 x 255 x 50 mm
B3-6625-70	665 x 255 x 70 mm
B3-6625-100	665 x 255 x 100 mm

Lid Code	Outside Dimension L x W x H
B3-6625	635 x 240 x 5 mm

Wide Size

Box Code	Outside Dimension L x W x H
B3-5333-50	533 x 337 x 50 mm
B3-5333-70	533 x 337 x 70 mm
B3-5333-100	533 x 337 x 100 mm

Lid Code	Outside Dimension L x W x H
B3-5333	503 x 322 x 5 mm

Fine Mesh Baskets

Full Size

Box Code	Outside Dimension L x W x H
B3-5425-50	540 x 255 x 50 mm
B3-5425-70	540 x 255 x 70 mm
B3-5425-100	540 x 255 x 100 mm

Lid Code	Outside Dimension L x W x H
B3-5425	510 x 240 x 5 mm

Din Size

Box Code	Outside Dimension L x W x H
B3-4825-50	485 x 255 x 50 mm
B3-4825-70	485 x 255 x 70 mm
B3-4825-100	485 x 255 x 100 mm

Lid Code	Outside Dimension L x W x H
B3-4825	455 x 240 x 5 mm

Quarter Size

Box Code	Outside Dimension L x W x H
B3-4025-50	405 x 255 x 50 mm
B3-4025-70	405 x 255 x 70 mm
B3-4025-100	405 x 255 x 100 mm

Lid Code	Outside Dimension L x W x H
B3-4025	375 x 240 x 5 mm

Ø 0,5 mm
□ 2 mm


Box Code	Outside Dimension L x W x H
B4-33-30	30 x 30 x 30 mm
B4-44-20	40 x 40 x 20 mm
B4-66-20	60 x 60 x 20 mm
B4-66-40	60 x 60 x 40 mm
B4-84-20	80 x 40 x 40 mm
B4-84-40	80 x 40 x 40 mm

Box Code	Outside Dimension L x W x H
B4-88-40	80 x 80 x 40 mm
B4-88-60	80 x 80 x 60 mm
B4-1212-30	120 x 120 x 30 mm
B4-1212-60	120 x 120 x 60 mm
B4-2424-60	240 x 240 x 60 mm

Scope Baskets



Perforated Sheet
Small Size Baskets



 Bottom With Lid
W/O Feet
Wire: $\varnothing$ 1 mm
Mesh: 4x4 mm

Code	Dimension L x W x H	Code	Dimension L x W x H
B2-2718-32	275 x 180 x 32 mm	B2-188-32	188 x 86 x 32 mm
B2-2718-42	275 x 180 x 42 mm	B2-144-28	146 x 46 x 28 mm
B2-247-30	240 x 76 x 30 mm	B2-84-28	86 x 46 x 28 mm
B2-1913-34	191 x 130 x 34 mm		

STU - SPRI Baskets

Wire Basket Single Scope	Wire Basket Double Scope	Perforated Basket Single Scope	Perforated Basket Double Scope
Code B5-298-55-S Dimension L x W x H 290 x 80 x 52 mm Code B5-468-55-S Dimension L x W x H 460 x 80 x 52 mm	Code B5-298-55-D Dimension L x W x H 290 x 80 x 52 mm Code B5-468-55-D Dimension L x W x H 460 x 80 x 52 mm	Code B5-278-44-S Dimension L x W x H 272 x 83 x 44 mm Code B5-458-44-S Dimension L x W x H 452 x 83 x 44 mm	Code B5-2710-44-D Dimension L x W x H 272 x 103 x 44 mm Code B5-4510-44-D Dimension L x W x H 452 x 103 x 44 mm

STU Baskets			
	Box Code	Outside Dimension L x W x H	Details
	B7-5728-135	575 x 280 x 135 mm	Upper Frame: $\varnothing$ 6 mm Body Wires: $\varnothing$ 2 & 3 mm Base Frame $\varnothing$ 5 mm
	5728-265	265 mm	
SPRI Baskets			
	Box Code	Outside Dimension L x W x H	Details
	B8-6030-40	600x300x40 mm	Main Frame: $\varnothing$ 6 mm
	B8-6030-100	600x300x100 mm	Body Wire: $\varnothing$ 4 mm
	B8-6030-150	600x300x150 mm	Base Frame $\varnothing$ 5 mm
	B8-6040-100	600x400x100 mm	Single Frame Stackable
	B8-6040-150	600x400x150 mm	
	B8-6040-200	600x400x200 mm	






Paper Filters		PTFE Teflon Filters	
Code A1-10-1 For: 1/1, 3/4, 1/2, Long, Wide Size: ø 185 mm	Code A1-11-1 For: Small Size Containers Size: ø 88 mm	Code A1-20-1 For: 1/1, 3/4, 1/2, Long, Wide Size: ø 190 mm Cycle: 2000 Cycles Reusable	Code A1-21-1 For: Small Size Containers Size: ø 88 mm Cycle: 2000 Cycles Reusable
Material 73g/m² Medical Grade Paper		Code A1-20-2 For: 1/1, 3/4, 1/2, Long, Wide Size: ø 190 mm Cycle: 1200 Cycles Reusable	Code A1-21-2 For: Small Size Containers Size: ø 88 mm Cycle: 1200 Cycles Reusable
Package 100 pcs			
Sterilization Steam			
Cycle Single Use			

Polymer Security Seals		Polymer Security Seals With Indicator	
Code A2-10-1	Code A2-20-1	Code A2-10-2	Code A2-20-2
Package 100 pcs	Package 100 pcs	Package 100 pcs	Package 100 pcs
Sterilization Steam	Sterilization Steam	Sterilization Steam	Sterilization Steam
Cycle Single Use	Cycle Single Use	Cycle Single Use	Cycle Single Use
Color Blue	Color Red	Color Blue	Color Red



Identification Labels				Filter Retainers	
Color	Code	Color	Code	Material 1,5 mm Aluminum	
	A5-10-		A5-10-	16 x 25 mm	
	1 A5-		5 A5-	Package 1 Piece	
	10-2		10-6	Sterilization Steam	
	A5-10-		A5-10-	Marking Maximum 24 Letters	
	3 A5-		7 A5-		
10-4 10-8 Paper Labels				Code A6-10-1 For: 1/1, 3/4, 1/2, Long, Wide	Code A6-10-2 For: Small Size Containers
	Code A3-10-1 with indicator		Code A3-10-2 with indicator	Material Stainless Steel	
				Package 1 piece	

Silicone Mats & Nets					
	Code	Size		Code	Size
	A7-	500 X 235 mm		A4-	540 x 252 mm
	5425	370 X 235 mm		5425	420 x 252 mm
	A7-	225 x 235 mm		A4-	254 x 253 mm
	4025	270 x 125 mm		4025	
A7-			A4-		
					
h: 16 mm			h: 16 mm		
Material			Material		
Package			Package		
Sterilization			Sterilization		
Color			Color		

Silicone Gasket					
	Code A11-10-1	Size 1/1 Full	Code A11-70-1	Size Mini	Material Medical Grade Sealing Silicone
	A11-20-13/4	Quarter	A11-50-1	Wide	Package 1 Piece
	A11-30-1	1/2 Half	A11-60-1	Dental	Shape D shape, Jointless
					Feature Place & Go, No Glue Required



Code T1-10-1
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece
Dimensions: 600 x 860 x 900 mm (W x L x H)



Code T1-10-2
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece
Dimensions: 550 x 860 x 900 mm (W x L x H)



Code T1-10-5
Material : Stainless Steel
Feet : 4 Adjustable Leveling Feet
Package : 1 Piece
Dimensions: 600 x 1400 x 850 mm (W x L x H)




Code T1-10-6
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece
Dimensions: 700 x 360 x 1350 mm (W x L x H)



Code T1-10-3
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece
Dimensions: 600 x 850 x 900 mm (W x L x H)



Code T1-10-4
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece
Dimensions: 550 x 860 x 900 mm (W x L x H)



Code T1-10-7
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece
Dimensions: 700 x 1080 x 1350 mm (W x L x H)



Code T1-10-8
Material : Stainless Steel
Feet : 4 wheels - 2 braked
Package : 1 Piece - including 3 baskets
Dimensions: 400 x 600 x 1700 (W x L x H)