

# Sterilization Container Systems

Secure Today  
Elevate Tomorrow







## Different Lids Colors

Our sterilization containers come in eight different colors (silver, yellow, red, blue, green, black, bronze and purple) ensuring to easily distinguish between different sets or types of instruments at a glance.




### Advantages of Multiple Colors

- **Enhanced Organization:** Different colors can be used to categorize and organize instruments based on their usage, department, or level of sterilization, reducing the risk of cross-contamination and improving efficiency..
- **Easy Identification:** Color-coding helps in quickly identifying the contents of each container, saving time and reducing errors in busy environments.
- **Improved Workflow:** With a clear visual system in place, staff can more easily manage inventory and ensure that the right instruments are always in the right place.

Customized Solutions & Tailored to Your Needs

Our sterilization containers are manufactured in 8 main sizes to accommodate various surgical needs: standard sizes (1/1 Full Size, 3/4 Quarter Size, 1/2 Half Size, Mini, Dental and 1/2 Dental Size) are ideal for common instrument sets, while the Wide Body size suits specialized instruments like Bookwalter retractors, and the Long Body size is designed for extended instruments such as robotic arms.

The sterilization containers are crafted in five primary heights to suit a variety of surgical instruments. Additionally, custom sizes or heights can be tailored to meet specific requirements.



Standard 8 Different Size Types  
For Standard & Reusable Models

| Type                                 | Packaging according to EN 868-8                                                                 |                    | Steam sterilization procedure according to EN 285 |            |                   |                 |                 |
|--------------------------------------|-------------------------------------------------------------------------------------------------|--------------------|---------------------------------------------------|------------|-------------------|-----------------|-----------------|
|                                      | Description                                                                                     | Symbolic Packaging | Gravity                                           | Pre-Vacuum | Fractional Vacuum | Fractional Flow | Steam Injection |
| Sterilization Container with filters | Textiles in container with filter in lid and box                                                |                    | +                                                 | +          | +                 | +               | +               |
|                                      |                                                                                                 |                    | O                                                 | O          | +                 | +               | +               |
|                                      | Textiles in container with filter in lid, non-perforated box                                    |                    | -                                                 | -          | +                 | +               | +               |
|                                      |                                                                                                 |                    | -                                                 | -          | +                 | +               | O               |
|                                      | Textiles in container with filter in box, non-perforated lid                                    |                    | O                                                 | +          | +                 | +               | +               |
|                                      |                                                                                                 |                    | -                                                 | O          | +                 | +               | +               |
|                                      | Instruments in wire basket, in sterilization filter in lid and box                              |                    | +3                                                | +          | +                 | +               | +               |
|                                      | Instruments in wire basket, in sterilization container filter in lid and non-perforated box     |                    | -                                                 | +2         | +2                | +2              | +2              |
|                                      | Instruments in wire basket, in sterilization container filter in box and non-perforated lid     |                    | O                                                 | +          | +                 | +               | +               |
|                                      |                                                                                                 |                    | -                                                 | +          | +                 | +               | +               |
| Sterilization Container with valves  | Textiles in sterilization container with valves in lid and box                                  |                    | -                                                 | O          | +                 | +               | O               |
|                                      |                                                                                                 |                    | -                                                 | -          | +                 | +               | O               |
|                                      | Textiles in sterilization container with valves in lid and no valves on the box                 |                    | -                                                 | -          | +                 | +               | -               |
|                                      |                                                                                                 |                    | -                                                 | -          | +                 | +               | -               |
|                                      | Instruments in wire basket in sterilization container with valve in lid and box                 |                    | -                                                 | +2         | +2                | +2              | +2              |
|                                      | Instruments in wire basket, in sterilization container with valve in lid and non-perforated box |                    | -                                                 | +2         | +2                | +2              | +2              |

\*To prevent errors during service, operation, or loading of the autoclaves, consult above table.  
\*Use only sterilization containers that guarantee proper evaporation or discharge of condensate.  
\*Keep in mind that gravitation-based procedures, which lack a dedicated drying process, do not ensure sufficient after-drying effects.



## STANDARD & REUSABLE MODELS

SECURE TODAY

ELEVATETOMORROW



## WHY STERILIZATION CONTAINERS?

- Consistent Processing
- Reduced Risk of Contamination
- Physical Instrument Protection
- Optimal Set Organization
- Reusability / Sustainability
- Reduced Waste
- Streamlined Process
- Easy Identification and Access
- Adherence to Standards
- Minimized Handling

# Today's Safety. Tomorrow's Efficiency



 Different Lids Colors

Danshealth.com

## Elevated Sterilization Solutions: For Your Operational Excellence

Discover our 2mm special alloy aluminum container systems that provide superior protection for your surgical instruments. Available in 8 different color palettes, they facilitate organization and quick identification. This solution supports your operational excellence by combining reliability, an eco-friendly approach, and long-term cost advantages.

### SAFETY LIDS



### TWO DIFFERENT FILTER TYPES

Paper Filter or PTFE Filter



### EASE WITH COLORS

EXCELLENCE IN WORKFLOW

 Different Lids Colors



### TWO BOTTOM TYPES

 Non Perforated

 Perforated



# Standard Model

1/1 Full Size

Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



Perforated



Code Outside Dimension Inside Dimension L x W x H

K1-10-110 597 x 275 x 95 mm 546 x 258 x 90 mm

K1-10-135 597 x 275 x 120 mm 546 x 258 x 115 mm

K1-10-150 597 x 275 x 136 mm 546 x 258 x 131 mm

K1-10-200 597 x 275 x 188 mm 546 x 258 x 183 mm

K1-10-260 597 x 275 x 235 mm 546 x 258 x 230 mm

Code Outside Dimension Inside Dimension L x W x H


KP1-10-110 597 x 275 x 95 mm 544 x 257 x 75 mm

KP1-10-135 597 x 275 x 120 mm 544 x 257 x 100 mm

KP1-10-150 597 x 275 x 136 mm 544 x 257 x 116 mm

KP1-10-200 597 x 275 x 188 mm 544 x 257 x 168 mm

KP1-10-260 597 x 275 x 235 mm 544 x 257 x 215 mm



Lids



Safety Lids



Colors Code

Grey K1-10-1

Yellow K1-10-2

Red K1-10-3

Blue K1-10-4

Green K1-10-5

Black K1-10-6

Purple K1-10-7

Bronze K1-10-8

Dimension 585 x 285 x 37 mm

Colors Code

Grey KS1-10-1

Yellow KS1-10-2

Red KS1-10-3

Blue KS1-10-4

Green KS1-10-5

Black KS1-10-6

Purple KS1-10-7


Bronze KS1-10-8

Dimension 586 x 289 x 36 mm

## Wire Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-5425-30  | 540 x 245 x 30 mm           |
| B1-5425-50  | 540 x 245 x 50 mm           |
| B1-5425-70  | 540 x 245 x 70 mm           |
| B1-5425-100 | 540 x 245 x 100 mm          |



Lid

Code B1-5425  
Outside Dimension L x W x H 510 x 230 x 5 mm

## Perforated Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-5425-30  | 540 x 251 x 40 mm           |
| B2-5425-50  | 540 x 251 x 60 mm           |
| B2-5425-70  | 540 x 251 x 80 mm           |
| B2-5425-100 | 540 x 251 x 110 mm          |



Lid

Code B2-5425  
Outside Dimension L x W x H 540 x 256 x 5 mm

## Hybrid Basket



Mesh: 4x4 mm Sides: ø 3 mm Wire: ø 1 mm



Box

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B3-5425-50  | 540 x 255 x 50 mm           |
| B3-5425-70  | 540 x 255 x 70 mm           |
| B3-5425-100 | 540 x 255 x 100 mm          |



Lid

Code B3-5425  
Outside Dimension L x W x H 510 x 240 x 5 mm

## Silicone Mat



Code A7-5425  
Outside Dimension L x W x H 500 x 235 mm

# Standard Model

1/1 Full Size Baskets

## Wire Basket



## Perforated Basket



## Hybrid Basket



# Standard Model

3/4Quarter Size

Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



Code Outside Dimension Inside Dimension L x W x H

K1-20-110 477 x 275 x 95 mm 423 x 258 x 90 mm

K1-20-135 477 x 275 x 120 mm 423 x 258 x 115 mm

K1-20-150 477 x 275 x 136 mm 423 x 258 x 131 mm

K1-20-200 477 x 275 x 188 mm 423 x 258 x 183 mm



Perforated




Code Outside Dimension Inside Dimension L x W x H

KP1-20-110 477 x 275 x 95 mm 423 x 258 x 75 mm

KP1-20-135 477 x 275 x 120 mm 423 x 258 x 100 mm

KP1-20-150 477 x 275 x 136 mm 423 x 258 x 116 mm

KP1-20-200 477 x 275 x 188 mm 423 x 258 x 168 mm



Lids



Colors Code

Grey K1-20-1  
Yellow K1-20-2  
Red K1-20-3  
Blue K1-20-4  
Green K1-20-5  
Black K1-20-6  
Purple K1-20-7  
Bronze K1-20-8

Dimension 465 x 280 x 37 mm



Safety Lids



Colors Code

Grey KS1-20-1  
Yellow KS1-20-2  
Red KS1-20-3  
Blue KS1-20-4  
Green KS1-20-5  
Black KS1-20-6  
Purple KS1-20-7  
Bronze KS1-20-8

Dimension 465 x 289 x 36 mm

## Wire Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box

Code Outside Dimension L x W x H  
B1-4025-30 405 x 245 x 30 mm  
B1-4025-50 405 x 245 x 50 mm  
B1-4025-70 405 x 245 x 70 mm  
B1-4025-100 405 x 245 x 100 mm




Lid

Code B1-4025  
Outside Dimension L x W x H 375 x 230 x 5 mm

## Perforated Basket



Mesh: 4x4 mm Wire: ø 1 mm



Box (with 10 mm feet)

Code Outside Dimension L x W x H  
B2-4025-30 403 x 251 x 40 mm  
B2-4025-50 403 x 251 x 60 mm  
B2-4025-70 403 x 251 x 80 mm  
B2-4025-100 403 x 251 x 110 mm




Lid

Code B2-4025  
Outside Dimension L x W x H 405 x 256 x 5 mm

## Hybrid Basket



Mesh: 4x4 mm Sides: ø 3 mm Wire: ø 1 mm



Box

Code Outside Dimension L x W x H  
B3-4025-50 405 x 255 x 50 mm  
B3-4025-70 405 x 255 x 70 mm  
B3-4025-100 405 x 255 x 100 mm



Lid

Code B3-4025  
Outside Dimension L x W x H 375 x 240 x 5 mm

## Silicone Mat




Code A7-4025  
Outside Dimension L x W x H 370 x 235 mm


# Standard Model

3/4 Quarter SizeBaskets


## Wire Basket



## Perforated Basket



## Hybrid Basket



## Standard Model

1/2 Half Size

Different Lids Colors



With Safety Lid



With Perforated Lid



Non Perforated



Code      Outside Dimension  
Inside Dimension  
L x W x H

K1-30-110    307 x 275 x 95 mm  
258 x 258 x 90 mm

K1-30-135    307 x 275 x 120 mm  
258 x 258 x 115 mm

K1-30-150    307 x 275 x 136 mm  
258 x 258 x 131 mm

K1-30-200    307 x 275 x 188 mm  
258 x 258 x 183 mm

K1-30-260    307 x 275 x 235 mm  
258 x 258 x 230 mm



Perforated



Code      Outside Dimension  
Inside Dimension  
L x W x H


KP1-30-110    307 x 275 x 95 mm  
258 x 258 x 75 mm

KP1-30-135    307 x 275 x 120 mm  
258 x 258 x 100 mm

KP1-30-150    307 x 275 x 136 mm  
258 x 258 x 116 mm

KP1-30-200    307 x 275 x 188 mm  
258 x 258 x 168 mm

KP1-30-260    307 x 275 x 235 mm  
258 x 258 x 215 mm



Lids



Colors      Code

Grey      K1-30-1

Yellow    K1-30-2

Red      K1-30-3

Blue     K1-30-4

Green    K1-30-5

Black    K1-30-6

Purple   K1-30-7

Bronze   K1-30-8

Dimension  
295 x 280 x 37 mm



Safety Lids



Colors      Code

Grey      KS1-30-1

Yellow    KS1-30-2

Red      KS1-30-3

Blue     KS1-30-4

Green    KS1-30-5

Black    KS1-30-6

Purple   KS1-30-7

Bronze   KS1-30-8

Dimension  
295 x 289 x 36 mm

### Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-2525-30  | 246 x 246 x 30 mm           |
| B1-2525-50  | 246 x 246 x 50 mm           |
| B1-2525-70  | 246 x 246 x 70 mm           |
| B1-2525-100 | 246 x 246 x 100 mm          |




Lid

| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| B1-2525 | 240 x 215 x 5 mm            |

### Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-2525-30  | 244 x 250 x 40 mm           |
| B2-2525-50  | 244 x 250 x 60 mm           |
| B2-2525-70  | 244 x 250 x 80 mm           |
| B2-2525-100 | 244 x 250 x 110 mm          |




Lid

| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| B2-2525 | 247 x 256 x 5 mm            |

### Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B3-2525-50  | 245 x 255 x 50 mm           |
| B3-2525-70  | 245 x 255 x 70 mm           |
| B3-2525-100 | 245 x 255 x 100 mm          |



Lid

| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| B3-2525 | 215 x 240 x 5 mm            |

### Silicone Mat



| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| A7-2525 | 225 x 235 mm                |

## Standard Model

1/2 Half Size Baskets

### Wire Basket



### Perforated Basket



### Hybrid Basket



# Standard Model

Long Size

Different Lids  
Colors



With Perforated Lid



With Color Options



Non  
Perforated



Code      Outside Dimension  
             Inside Dimension  
             L x W x H

K1-40-150    770 x 275 x 136 mm  
                 715 x 258 x 131 mm

K1-40-200    770 x 275 x 188 mm  
                 715 x 258 x 183 mm



Perforated



Code      Outside Dimension  
             Inside Dimension  
             L x W x H

KP1-40-150    770 x 275 x 136 mm  
                 715 x 258 x 116 mm

KP1-40-200    770 x 275 x 188 mm  
                 715 x 258 x 168 mm

Lids



Colors      Code

Grey      K1-40-1

Yellow    K1-40-2

Red      K1-40-3

Blue      K1-40-4

Green    K1-40-5

Black    K1-40-6

Purple    K1-40-7

Bronze    K1-40-8

Dimension  
758 x 280 x 37 mm

## Wire Basket



Box

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B1-6625-50  | 665 x 245 x 50 mm           |
| B1-6625-70  | 665 x 245 x 70 mm           |
| B1-6625-100 | 665 x 245 x 100 mm          |

Mesh: 4x4 mm  
Wire: ø 1 mm

Lid

| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| B1-6625 | 635 x 230 x 5 mm            |

## Perforated Basket



Box (with 10 mm feet)

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B2-6625-30  | 650 x 251 x 40 mm           |
| B2-6625-50  | 650 x 251 x 60 mm           |
| B2-6625-70  | 650 x 251 x 80 mm           |
| B2-6625-100 | 650 x 251 x 110 mm          |

Mesh: 4x4 mm  
Wire: ø 1 mm

Lid

| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| B2-6625 | 652 x 256 x 5 mm            |

## Hybrid Basket



Box

| Code        | Outside Dimension L x W x H |
|-------------|-----------------------------|
| B3-6625-50  | 665 x 255 x 50 mm           |
| B3-6625-70  | 665 x 255 x 70 mm           |
| B3-6625-100 | 665 x 255 x 100 mm          |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm

Lid

| Code    | Outside Dimension L x W x H |
|---------|-----------------------------|
| B3-6625 | 635 x 240 x 5 mm            |


# Standard Model

Long Size Baskets

## Wire Basket



## Perforated Basket



## Hybrid Basket



# Standard Model

Wide Size

 Different Lids  
Colors



With Perforated Lid



With Color Options



Non  
Perforated



Code      Outside Dimension  
             Inside Dimension  
             L x W x H

K1-50-110    597 x 365 x 105 mm  
                 546 x 348 x 100 mm

K1-50-170    597 x 365 x 165 mm  
                 546 x 348 x 160 mm




Perforated



Code      Outside Dimension  
             Inside Dimension  
             L x W x H

KP1-50-110    597 x 365 x 105 mm  
                 546 x 348 x 85 mm

KP1-50-170    597 x 365 x 165 mm  
                 546 x 348 x 145 mm

 Lids

Colors      Code

Grey      K1-50-1

Yellow    K1-50-2

Red      K1-50-3

Blue      K1-50-4

Green    K1-50-5

Black    K1-50-6

Purple    K1-50-7

Bronze    K1-50-8

Dimension  
585 x 370 x 37 mm



## Wire Basket



 Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B1-5333-50  | 533 x 337 x 50 mm              |
| B1-5333-70  | 533 x 337 x 70 mm              |
| B1-5333-100 | 533 x 337 x 100 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm

 Lid

Code  
B1-5333  
Outside Dimension  
L x W x H  
503 x 322 x 5 mm

## Perforated Basket



 Box (with 10 mm feet)

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B2-5333-30  | 535 x 340 x 40 mm              |
| B2-5333-50  | 535 x 340 x 60 mm              |
| B2-5333-70  | 535 x 340 x 80 mm              |
| B2-5333-100 | 535 x 340 x 110 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm

 Lid

Code  
B2-5333  
Outside Dimension  
L x W x H  
537 x 343 x 5 mm

## Hybrid Basket



 Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B3-5333-50  | 533 x 337 x 50 mm              |
| B3-5333-70  | 533 x 337 x 70 mm              |
| B3-5333-100 | 533 x 337 x 100 mm             |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm

 Lid

Code  
B3-5333  
Outside Dimension  
L x W x H  
503 x 322 x 5 mm

# Standard Model

Wide Size Baskets

## Wire Basket



## Perforated Basket



## Hybrid Basket



# Standard Model

## Flat Size


 Different Lids  
Colors



With Perforated Lid



With Color Options



Non  
Perforated

Code      Outside Dimension  
             Inside Dimension  
             L x W x H

|          |                                        |
|----------|----------------------------------------|
| K1-90-45 | 297 x 275 x 55 mm<br>258 x 258 x 50 mm |
| K1-90-75 | 297 x 275 x 75 mm<br>258 x 258 x 70 mm |




Lids


| Colors | Code    |
|--------|---------|
| Grey   | K1-90-1 |
| Yellow | K1-90-2 |
| Red    | K1-90-3 |
| Blue   | K1-90-4 |
| Green  | K1-90-5 |
| Black  | K1-90-6 |
| Purple | K1-90-7 |
| Bronze | K1-90-8 |

Dimension  
295 x 280 x 37 mm

## Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code       | Outside Dimension<br>L x W x H |
|------------|--------------------------------|
| B1-2525-30 | 246 x 246 x 30 mm              |
| B1-2525-50 | 246 x 246 x 50 mm              |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2525 | 240 x 215 x 5 mm               |

## Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code       | Outside Dimension<br>L x W x H |
|------------|--------------------------------|
| B2-2525-30 | 244 x 250 x 40 mm              |
| B2-2525-50 | 244 x 250 x 60 mm              |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2525 | 247 x 256 x 5 mm               |


## Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2525 | 225 x 235 mm                   |

# Standard Model

## Flat Size Baskets



Wire Basket

## Perforated Basket



# Standard Model

## MiniSize

 Different Lids  
Colors



With Perforated Lid



With Color Options



Non  
Perforated

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-70-40  | 305 x 135 x 45 mm<br>272 x 125 x 33 mm             |
| K1-70-70  | 305 x 135 x 60 mm<br>272 x 125 x 48 mm             |
| K1-70-100 | 305 x 135 x 90 mm<br>272 x 125 x 78 mm             |
| K1-70-130 | 305 x 135 x 120 mm<br>272 x 125 x 108 mm           |



Lids

| Colors | Code    |
|--------|---------|
| Grey   | K1-70-1 |
| Yellow | K1-70-2 |
| Red    | K1-70-3 |
| Blue   | K1-70-4 |
| Green  | K1-70-5 |
| Black  | K1-70-6 |
| Purple | K1-70-7 |
| Bronze | K1-70-8 |

Dimension  
303 x 136 x 29 mm

## Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2712-40  | 268 x 117 x 40 mm              |
| B1-2712-100 | 268 x 117 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2712 | 250 x 120 x 5 mm               |

## Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2712-40  | 260 x 120 x 50 mm              |
| B2-2712-100 | 260 x 120 x 110 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2712 | 262 x 122 x 5 mm               |

## Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2712 | 270 x 125 mm                   |


# Standard Model

## MiniSizeBaskets

### Wire Basket



### Perforated Basket



# Standard Model

Dental Size

 Different Lids  
Colors



With Perforated Lid




With Color Options



Non  
Perforated

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-60-40  | 310 x 189 x 45 mm<br>277 x 180 x 33 mm             |
| K1-60-70  | 310 x 189 x 60 mm<br>277 x 180 x 48 mm             |
| K1-60-100 | 310 x 189 x 90 mm<br>277 x 180 x 78 mm             |
| K1-60-130 | 310 x 189 x 120 mm<br>277 x 180 x 108 mm           |



Lids

| Colors | Code    |
|--------|---------|
| Grey   | K1-60-1 |
| Yellow | K1-60-2 |
| Red    | K1-60-3 |
| Blue   | K1-60-4 |
| Green  | K1-60-5 |
| Black  | K1-60-6 |
| Purple | K1-60-7 |
| Bronze | K1-60-8 |

Dimension  
308 x 190 x 29 mm

## Wire Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2717-40  | 270 x 170 x 40 mm              |
| B1-2717-100 | 270 x 170 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2717 | 252 x 172 x 5 mm               |

## Perforated Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2717-40  | 273 x 172 x 50 mm              |
| B2-2717-100 | 273 x 172 x 110 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2717 | 275 x 174 x 5 mm               |


# Standard Model

Dental Size Baskets

## Wire Basket



## Perforated Basket



# Standard Model

1/2 Dental Size

 Different Lids  
Colors



With Perforated Lid



With Color Options



Non  
Perforated

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-80-40  | 207 x 144 x 45 mm<br>174 x 135 x 33 mm             |
| K1-80-70  | 207 x 144 x 60 mm<br>174 x 135 x 48 mm             |
| K1-80-100 | 207 x 144 x 90 mm<br>174 x 135 x 78 mm             |
| K1-80-130 | 207 x 144 x 120 mm<br>174 x 135 x 108 mm           |

 Lids

| Colors | Code    |
|--------|---------|
| Grey   | K1-80-1 |
| Yellow | K1-80-2 |
| Red    | K1-80-3 |
| Blue   | K1-80-4 |
| Green  | K1-80-5 |
| Black  | K1-80-6 |
| Purple | K1-80-7 |
| Bronze | K1-80-8 |

Dimension  
205 x 145 x 29 mm

## Wire Basket



Mesh:  $\varnothing$  4x4 mm  
Wire:  $\varnothing$  1 mm



Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-1714-40  | 165 x 130 x 35 mm              |
| B1-1714-100 | 165 x 130 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-1714 | 147 x 132 x 5 mm               |

## Perforated Basket



Mesh:  $\varnothing$  4x4 mm  
Wire:  $\varnothing$  1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-1714-40  | 165 x 130 x 45 mm              |
| B2-1714-100 | 165 x 130 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-1714 | 167 x 130 x 5 mm               |

# Standard Model

1/2 Dental Size Baskets

## Wire Basket



## Perforated Basket



# Clean. Safe. Sustainable.

2 Different Filtering System  
Mechanical Valve System & Barrier System




Different Lids Colors

Danhealth.com

## BARRIER MODEL

The Barrier Model is an innovative sterilization container designed to ensure maximum safety in CSSD workflows. Its labyrinth structure prevents microorganisms from entering, maintaining sterility for extended periods. This design offers reliable, practical, and scientifically proven protection.


BARRIER FILTERS



EASY-TO-USE LOCK / UNLOCK



STEAM FILTER CROSS-SECTIONAL VIEW



EASE WITH COLORS  
EXCELLENCE IN WORKFLOW

Different Lids Colors



### Scientific Foundation

The principle behind the Barrier mechanism traces back to Louis Pasteur's groundbreaking "Swan Neck Flask" experiment in 1862, which challenged the idea of spontaneous generation. Pasteur designed a flask with an S-shaped neck that allowed air to pass freely but trapped dust and microorganisms in its curves. When the broth inside was boiled and left untouched, no microbial growth occurred, clearly proving that life could not emerge without contamination from existing microorganisms. However, once the flask was tilted and the trapped particles reached the broth, microbial growth appeared almost immediately.

The Barrier Model adopts the same principle through its labyrinth-like system. By forcing air to travel through a complex path, it prevents microorganisms and viruses from penetrating the container. This ensures that sterility is maintained, while still allowing necessary airflow. The result is a reliable and scientifically validated protective system, directly inspired by one of the most important experiments in microbiology.



REDUCE. REUSE. RENEW.

Seize Cost Control and Sustainability: Ditch the waste expense and embrace the smart, permanent solution.



VALVE MODEL

Valve Model sterilization containers are designed with an advanced mechanical valve system that removes the need for disposable filters while delivering stable and reliable performance throughout extended use.

Each valve is precisely engineered with integrated gaskets that ensure smooth operation in response to pressure variations during the sterilization process. The secure closure, supported by high-grade silicone gaskets, provides excellent sealing integrity and long-term durability.


Combining innovative design with practicality, Dans Valve Model containers offer a sustainable and cost-efficient solution for safe and effective sterilization in healthcare environments worldwide.

EASE WITH COLORS EXCELLENCE IN WORKFLOW

Different Lids Colors



RUNNING PRINCIPLE



At the beginning of the sterilization cycle, the vacuum phase causes the valves to fully open, allowing the evacuation of dry air from the container.

During the sterilization cycle, the pressure-activated valves automatically open under positive pressure, enabling controlled steam penetration for effective sterilization.

Once the process is complete, and no external forces such as vacuum or pressure are present, the entire system remains sealed, preserving the sterility of its contents.

# Reusable Model

1/1 Full Size

 Different Lids  
Colors



With Reusable Mechanism



With Color Options



Non  
Perforated


| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-10-110 | 597 x 275 x 95 mm<br>546 x 258 x 90 mm             |
| K1-10-135 | 597 x 275 x 120 mm<br>546 x 258 x 115 mm           |
| K1-10-150 | 597 x 275 x 136 mm<br>546 x 258 x 131 mm           |
| K1-10-200 | 597 x 275 x 188 mm<br>546 x 258 x 183 mm           |
| K1-10-260 | 597 x 275 x 235 mm<br>546 x 258 x 230 mm           |



Valve Lids

| Colors | Code    |
|--------|---------|
| Grey   | K2-10-1 |
| Yellow | K2-10-2 |
| Red    | K2-10-3 |
| Blue   | K2-10-4 |
| Green  | K2-10-5 |
| Black  | K2-10-6 |
| Purple | K2-10-7 |
| Bronze | K2-10-8 |

Dimension  
585 x 280 x 37 mm




Barrier Lids

| Colors | Code    |
|--------|---------|
| Grey   | K3-10-1 |
| Yellow | K3-10-2 |
| Red    | K3-10-3 |
| Blue   | K3-10-4 |
| Green  | K3-10-5 |
| Black  | K3-10-6 |
| Purple | K3-10-7 |
| Bronze | K3-10-8 |

Dimension  
585 x 280 x 37 mm

## Wire Basket



Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B1-5425-30  | 540 x 245 x 30 mm              |
| B1-5425-50  | 540 x 245 x 50 mm              |
| B1-5425-70  | 540 x 245 x 70 mm              |
| B1-5425-100 | 540 x 245 x 100 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-5425 | 510 x 230 x 5 mm               |

## Perforated Basket



Box (with 10 mm feet)

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B2-5425-30  | 540 x 251 x 40 mm              |
| B2-5425-50  | 540 x 251 x 60 mm              |
| B2-5425-70  | 540 x 251 x 80 mm              |
| B2-5425-100 | 540 x 251 x 110 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-5425 | 540 x 256 x 5 mm               |

## Hybrid Basket



Box

| Code        | Outside Dimension L x<br>W x H |
|-------------|--------------------------------|
| B3-5425-50  | 540 x 255 x 50 mm              |
| B3-5425-70  | 540 x 255 x 70 mm              |
| B3-5425-100 | 540 x 255 x 100 mm             |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-5425 | 510 x 240 x 5 mm               |

## Silicone Mat




| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-5425 | 500 x 235 mm                   |


# Reusable Model

1/1 Full Size Baskets


## Wire Basket



## Perforated Basket



## Hybrid Basket



# Reusable Model

3/4 Quarter Size

 Different Lids  
Colors



With Reusable Mechanism



With Color Options



| Code      | Outside Dimension                        |
|-----------|------------------------------------------|
|           | Inside Dimension<br>L x W x H            |
| K1-20-110 | 477 x 275 x 95 mm<br>423 x 258 x 90 mm   |
| K1-20-135 | 477 x 275 x 120 mm<br>423 x 258 x 115 mm |
| K1-20-150 | 477 x 275 x 136 mm<br>423 x 258 x 131 mm |
| K1-20-200 | 477 x 275 x 188 mm<br>423 x 258 x 183 mm |



| Colors | Code    |
|--------|---------|
| Grey   | K2-20-1 |
| Yellow | K2-20-2 |
| Red    | K2-20-3 |
| Blue   | K2-20-4 |
| Green  | K2-20-5 |
| Black  | K2-20-6 |
| Purple | K2-20-7 |
| Bronze | K2-20-8 |


Dimension  
465 x 280 x 37 mm



| Colors | Code    |
|--------|---------|
| Grey   | K3-20-1 |
| Yellow | K3-20-2 |
| Red    | K3-20-3 |
| Blue   | K3-20-4 |
| Green  | K3-20-5 |
| Black  | K3-20-6 |
| Purple | K3-20-7 |
| Bronze | K3-20-8 |

Dimension  
465 x 280 x 37 mm

## Wire Basket



| Code        | Outside Dimension L x<br>W x H |
|-------------|--------------------------------|
| B1-4025-30  | 405 x 245 x 30 mm              |
| B1-4025-50  | 405 x 245 x 50 mm              |
| B1-4025-70  | 405 x 245 x 70 mm              |
| B1-4025-100 | 405 x 245 x 100 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-4025 | 375 x 230 x 5 mm               |

## Perforated Basket




| Code        | Outside Dimension L x<br>W x H |
|-------------|--------------------------------|
| B2-4025-30  | 403 x 251 x 40 mm              |
| B2-4025-50  | 403 x 251 x 60 mm              |
| B2-4025-70  | 403 x 251 x 80 mm              |
| B2-4025-100 | 403 x 251 x 110 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-4025 | 405 x 256 x 5 mm               |

## Hybrid Basket



| Code        | Outside Dimension L x<br>W x H |
|-------------|--------------------------------|
| B3-4025-50  | 405 x 255 x 50 mm              |
| B3-4025-70  | 405 x 255 x 70 mm              |
| B3-4025-100 | 405 x 255 x 100 mm             |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-4025 | 375 x 240 x 5 mm               |

## Silicone Mat




| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-4025 | 370 x 235 mm                   |


# Reusable Model

3/4 Quarter Size Baskets


## Wire Basket



## Perforated Basket



## Hybrid Basket



# Reusable Model

1/2 Half Size

Different Lids Colors



With Reusable Mechanism



With Color Options



Non Perforated

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-30-110 | 307 x 275 x 95 mm<br>258 x 258 x 90 mm             |
| K1-30-135 | 307 x 275 x 120 mm<br>258 x 258 x 115 mm           |
| K1-30-150 | 307 x 275 x 136 mm<br>258 x 258 x 131 mm           |
| K1-30-200 | 307 x 275 x 188 mm<br>258 x 258 x 183 mm           |
| K1-30-260 | 307 x 275 x 235 mm<br>258 x 258 x 230 mm           |



Valve Lids

| Colors | Code    |
|--------|---------|
| Grey   | K2-30-1 |
| Yellow | K2-30-2 |
| Red    | K2-30-3 |
| Blue   | K2-30-4 |
| Green  | K2-30-5 |
| Black  | K2-30-6 |
| Purple | K2-30-7 |
| Bronze | K2-30-8 |

Dimension  
295 x 280 x 37 mm



Barrier Lids

| Colors | Code    |
|--------|---------|
| Grey   | K3-30-1 |
| Yellow | K3-30-2 |
| Red    | K3-30-3 |
| Blue   | K3-30-4 |
| Green  | K3-30-5 |
| Black  | K3-30-6 |
| Purple | K3-30-7 |
| Bronze | K3-30-8 |

Dimension  
295 x 280 x 37 mm

## Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B1-2525-30  | 246 x 246 x 30 mm              |
| B1-2525-50  | 246 x 246 x 50 mm              |
| B1-2525-70  | 246 x 246 x 70 mm              |
| B1-2525-100 | 246 x 246 x 100 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2525 | 240 x 215 x 5 mm               |

## Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B2-2525-30  | 244 x 250 x 40 mm              |
| B2-2525-50  | 244 x 250 x 60 mm              |
| B2-2525-70  | 244 x 250 x 80 mm              |
| B2-2525-100 | 244 x 250 x 110 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2525 | 247 x 256 x 5 mm               |

## Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B3-2525-50  | 245 x 255 x 50 mm              |
| B3-2525-70  | 245 x 255 x 70 mm              |
| B3-2525-100 | 245 x 255 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-2525 | 215 x 240 x 5 mm               |

## Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2525 | 225 x 235 mm                   |

# Reusable Model

1/2 Half Size Baskets

## Wire Basket



## Perforated Basket



## Hybrid Basket



# Reusable Model

Long Size

 Different Lids  
Colors



With Reusable Mechanism



With Color Options



| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-40-150 | 770 x 275 x 136 mm<br>715 x 258 x 131 mm           |
| K1-40-200 | 770 x 275 x 188 mm<br>715 x 258 x 183 mm           |




| Colors    | Code                      |
|-----------|---------------------------|
| Grey      | K2-40-                    |
| Yellow    | 1 K2-                     |
| Red       | 40-2                      |
| Blue      | K2-40-                    |
| Green     | 3 K2-                     |
| Black     | 40-4                      |
| Purple    | K2-40-                    |
| Bronze    | 5 K2-                     |
| Dimension | 40-6<br>758 x 280 x 37 mm |
|           | 7 K2-                     |
|           | 40-8                      |



| Colors    | Code                        |
|-----------|-----------------------------|
| Grey      | K3-40-                      |
| Yellow    | 1 K3-                       |
| Red       | 40-2                        |
| Blue      | K3-40-                      |
| Green     | 3 K3-                       |
| Black     | 40-4                        |
| Purple    | K3-40-                      |
| Bronze    | 5                           |
| Dimension | K3-40-<br>758 x 280 x 37 mm |
|           | K3-40-                      |
|           | 7                           |
|           | K3-40-                      |
|           | 8                           |

## Wire Basket



| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B1-6625-50  | 665 x 245 x 50 mm              |
| B1-6625-70  | 665 x 245 x 70 mm              |
| B1-6625-100 | 665 x 245 x 100 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-6625 | 635 x 230 x 5 mm               |

## Perforated Basket



| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B2-6625-30  | 650 x 251 x 40 mm              |
| B2-6625-50  | 650 x 251 x 60 mm              |
| B2-6625-70  | 650 x 251 x 80 mm              |
| B2-6625-100 | 650 x 251 x 110 mm             |

Mesh: 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-6625 | 652 x 256 x 5 mm               |

## Hybrid Basket



| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B3-6625-50  | 665 x 255 x 50 mm              |
| B3-6625-70  | 665 x 255 x 70 mm              |
| B3-6625-100 | 665 x 255 x 100 mm             |

Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-6625 | 635 x 240 x 5 mm               |

# Reusable Model

Long Size Baskets

## Wire Basket



## Perforated Basket



## Hybrid Basket



Reusable Model

Wide Size

Different Lids  
Colors



With Reusable Mechanism



With Color Options



Non  
Perforated

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-50-110 | 597 x 365 x 105 mm<br>546 x 348 x 100 mm           |
| K1-50-170 | 597 x 365 x 165 mm<br>546 x 348 x 160 mm           |



Valve Lids

| Colors | Code    |
|--------|---------|
| Grey   | K2-50-1 |
| Yellow | K2-50-2 |
| Red    | K2-50-3 |
| Blue   | K2-50-4 |
| Green  | K2-50-5 |
| Black  | K2-50-6 |
| Purple | K2-50-7 |
| Bronze | K2-50-8 |

Dimension  
585 x 370 x 37mm



Barrier Lids


| Colors | Code    |
|--------|---------|
| Grey   | K3-50-1 |
| Yellow | K3-50-2 |
| Red    | K3-50-3 |
| Blue   | K3-50-4 |
| Green  | K3-50-5 |
| Black  | K3-50-6 |
| Purple | K3-50-7 |
| Bronze | K3-50-8 |

Dimension  
585 x 370 x 37mm

Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B1-5333-50  | 533 x 337 x 50 mm              |
| B1-5333-70  | 533 x 337 x 70 mm              |
| B1-5333-100 | 533 x 337 x 100 mm             |



Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-5333 | 503 x 322 x 5 mm               |

Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)

| Code        | Outside Dimension L x<br>W x H |
|-------------|--------------------------------|
| B2-5333-30  | 535 x 340 x 40 mm              |
| B2-5333-50  | 535 x 340 x 60 mm              |
| B2-5333-70  | 535 x 340 x 80 mm              |
| B2-5333-100 | 535 x 340 x 110 mm             |




Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-5333 | 537 x 343 x 5 mm               |

Hybrid Basket



Mesh: 4x4 mm  
Sides: ø 3 mm  
Wire: ø 1 mm



Box

| Code        | Outside Dimension L<br>x W x H |
|-------------|--------------------------------|
| B3-5333-50  | 533 x 337 x 50 mm              |
| B3-5333-70  | 533 x 337 x 70 mm              |
| B3-5333-100 | 533 x 337 x 100 mm             |




Lid


| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B3-5333 | 503 x 322 x 5 mm               |

Reusable Model

Wide Size Baskets




Wire Basket



Perforated Basket



Hybrid Basket



# Reusable Model

Flat Size

 Different Lids  
Colors



With Reusable Mechanism



With Color Options




| Code     | Outside Dimension<br>Inside Dimension<br>L x W x H |
|----------|----------------------------------------------------|
| K1-90-75 | 297 x 275 x 75 mm<br>258 x 258 x 70 mm             |




| Colors | Code    |
|--------|---------|
| Grey   | K3-90-1 |
| Yellow | K3-90-2 |
| Red    | K3-90-3 |
| Blue   | K3-90-4 |
| Green  | K3-90-5 |
| Black  | K3-90-6 |
| Purple | K3-90-7 |
| Bronze | K3-90-8 |

Dimension  
295 x 280 x 37 mm

## Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



| Code       | Outside Dimension<br>L x W x H |
|------------|--------------------------------|
| B1-2525-30 | 246 x 246 x 30 mm              |
| B1-2525-50 | 246 x 246 x 50 mm              |




| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2525 | 240 x 215 x 5 mm               |

## Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



| Code       | Outside Dimension<br>L x W x H |
|------------|--------------------------------|
| B2-2525-30 | 244 x 250 x 40 mm              |
| B2-2525-50 | 244 x 250 x 60 mm              |



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2525 | 247 x 256 x 5 mm               |

## Silicone Mat




| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2525 | 225 x 235 mm                   |


# Reusable Model

Flat Size Baskets

## Wire Basket



## Perforated Basket



# Reusable Model

## Mini Size

 Different Lids  
Colors



With Reusable Mechanism



With Color Options



Non  
Perforated

| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-70-70  | 305 x 135 x 60 mm<br>272 x 125 x 48 mm             |
| K1-70-100 | 305 x 135 x 90 mm<br>272 x 125 x 78 mm             |
| K1-70-130 | 305 x 135 x 120 mm<br>272 x 125 x 108 mm           |



Lids

| Colors | Code    |
|--------|---------|
| Grey   | K3-70-1 |
| Yellow | K3-70-2 |
| Red    | K3-70-3 |
| Blue   | K3-70-4 |
| Green  | K3-70-5 |
| Black  | K3-70-6 |
| Purple | K3-70-7 |
| Bronze | K3-70-8 |

Dimension  
303 x 136 x 29 mm

## Wire Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box



Lid

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2712-40  | 268 x 117 x 40 mm              |
| B1-2712-100 | 268 x 117 x 100 mm             |

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2712 | 250 x 120 x 5 mm               |

## Perforated Basket



Mesh: 4x4 mm  
Wire: ø 1 mm



Box (with 10 mm feet)



Lid

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2712-40  | 260 x 120 x 50 mm              |
| B2-2712-100 | 260 x 120 x 110 mm             |

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2712 | 262 x 122 x 5 mm               |

## Silicone Mat



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| A7-2712 | 270 x 125 mm                   |


# Reusable Model

## Mini Size Baskets

### Wire Basket



### Perforated Basket



# Reusable Model

Dental Size

 Different Lids  
Colors



With Reusable Mechanism



With Color Options



| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-60-70  | 310 x 189 x 60 mm<br>277 x 180 x 48 mm             |
| K1-60-100 | 310 x 189 x 90 mm<br>277 x 180 x 78 mm             |
| K1-60-130 | 310 x 189 x 120 mm<br>277 x 180 x 108 mm           |



| Colors | Code    |
|--------|---------|
| Grey   | K3-60-1 |
| Yellow | K3-60-2 |
| Red    | K3-60-3 |
| Blue   | K3-60-4 |
| Green  | K3-60-5 |
| Black  | K3-60-6 |
| Purple | K3-60-7 |
| Bronze | K3-60-8 |

Dimension  
308 x 190 x 29 mm

## Wire Basket




| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2717-40  | 270 x 170 x 40 mm              |
| B1-2717-100 | 270 x 170 x 100 mm             |

Mesh: ø 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-2717 | 252 x 172 x 5 mm               |

## Perforated Basket



| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2717-40  | 273 x 172 x 50 mm              |
| B2-2717-100 | 273 x 172 x 110 mm             |

Mesh: ø 4x4 mm  
Wire: ø 1 mm



| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-2717 | 275 x 174 x 5 mm               |

# Reusable Model

Dental Size Baskets

## Wire Basket



## Perforated Basket



# Reusable Model

1/2 Dental Size

Different Lids Colors



With Reusable Mechanism



With Color Options



Non Perforated


| Code      | Outside Dimension<br>Inside Dimension<br>L x W x H |
|-----------|----------------------------------------------------|
| K1-80-70  | 207 x 144 x 60 mm<br>174 x 135 x 48 mm             |
| K1-80-100 | 207 x 144 x 90 mm<br>174 x 135 x 78 mm             |
| K1-80-130 | 207 x 144 x 120 mm<br>174 x 135 x 108 mm           |

Lids

| Colors | Code    |
|--------|---------|
| Grey   | K3-80-1 |
| Yellow | K3-80-2 |
| Red    | K3-80-3 |
| Blue   | K3-80-4 |
| Green  | K3-80-5 |
| Black  | K3-80-6 |
| Purple | K3-80-7 |
| Bronze | K3-80-8 |

Dimension  
205 x 145 x 29 mm

## Wire Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm

Box

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-1714-40  | 165 x 130 x 35 mm              |
| B1-1714-100 | 165 x 130 x 100 mm             |

Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B1-1714 | 147 x 132 x 5 mm               |



## Perforated Basket

## Perforated Basket



Mesh: ø 4x4 mm  
Wire: ø 1 mm

Box (with 10 mm feet)

| Code        | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-1714-40  | 165 x 130 x 45 mm              |
| B2-1714-100 | 165 x 130 x 110 mm             |

Lid

| Code    | Outside Dimension<br>L x W x H |
|---------|--------------------------------|
| B2-1714 | 167 x 130 x 5 mm               |




Organize.  
Sterilize.

Excel.


Holding the  
Standard




Three Different Types of Basket  
Models by Construction/Material



Perforated  
Sheet  
Ø 1 mm  
□ 4 mm



Wire Mesh  
Ø 1 mm  
□ 4 mm



Hybrith  
Ø 1 mm  
□ 4 mm  
Side Mesh  
Ø 3 mm



Half Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2525-30  | 246 x 246 x 30 mm              |
| B1-2525-50  | 246 x 246 x 50 mm              |
| B1-2525-70  | 246 x 246 x 70 mm              |
| B1-2525-100 | 246 x 246 x 100 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-2525  | 240 x 215 x 5 mm               |

Long Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-6625-50  | 665 x 245 x 50 mm              |
| B1-6625-70  | 665 x 245 x 70 mm              |
| B1-6625-100 | 665 x 245 x 100 mm             |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-6625  | 635 x 230 x 5 mm               |

Wide Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-5333-50  | 533 x 337 x 50 mm              |
| B1-5333-70  | 533 x 337 x 70 mm              |
| B1-5333-100 | 533 x 337 x 100 mm             |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-5333  | 503 x 322 x 5 mm               |

Full Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-5425-30  | 540 x 245 x 30 mm              |
| B1-5425-50  | 540 x 245 x 50 mm              |
| B1-5425-70  | 540 x 245 x 70 mm              |
| B1-5425-100 | 540 x 245 x 100 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-5425  | 510 x 230 x 5 mm               |

Din Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-4825-30  | 480 x 245 x 30 mm              |
| B1-4825-50  | 480 x 245 x 50 mm              |
| B1-4825-70  | 480 x 245 x 70 mm              |
| B1-4825-100 | 480 x 245 x 100 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-4825  | 450 x 230 x 5 mm               |

Quarter Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-4025-30  | 405 x 245 x 30 mm              |
| B1-4025-50  | 405 x 245 x 50 mm              |
| B1-4025-70  | 405 x 245 x 70 mm              |
| B1-4025-100 | 405 x 245 x 100 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-4025  | 375 x 230 x 5 mm               |

Dental Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2717-40  | 270 x 170 x 40 mm              |
| B1-2717-100 | 270 x 170 x 100 mm             |
| -           | -                              |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-2717  | 252 x 172 x 5 mm               |

Mini Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-2712-40  | 268 x 117 x 40 mm              |
| B1-2712-100 | 268 x 117 x 100 mm             |
| -           | -                              |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-2712  | 250 x 120 x 5 mm               |


Dental 1/2



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B1-1714-40  | 165 x 130 x 35 mm              |
| B1-1714-100 | 165 x 130 x 100 mm             |
| -           | -                              |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B1-1714  | 147 x 132 x 5 mm               |



Half Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2525-30  | 244 x 250 x 40 mm              |
| B2-2525-50  | 244 x 250 x 60 mm              |
| B2-2525-70  | 244 x 250 x 80 mm              |
| B2-2525-100 | 244 x 250 x 110 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-2525  | 247 x 256 x 5 mm               |

Long Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-6625-30  | 650 x 251 x 40 mm              |
| B2-6625-50  | 650 x 251 x 60 mm              |
| B2-6625-70  | 650 x 251 x 80 mm              |
| B2-6625-100 | 650 x 251 x 110 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-6625  | 652 x 256 x 5 mm               |

Wide Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-5333-30  | 535 x 340 x 40 mm              |
| B2-5333-50  | 535 x 340 x 60 mm              |
| B2-5333-70  | 535 x 340 x 80 mm              |
| B2-5333-100 | 535 x 340 x 110 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-5333  | 537 x 343 x 5 mm               |

Full Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-5425-30  | 540 x 251 x 40 mm              |
| B2-5425-50  | 540 x 251 x 60 mm              |
| B2-5425-70  | 540 x 251 x 80 mm              |
| B2-5425-100 | 540 x 251 x 110 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-5425  | 540 x 256 x 5 mm               |

Din Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-4825-30  | 483 x 251 x 40 mm              |
| B2-4825-50  | 483 x 251 x 60 mm              |
| B2-4825-70  | 483 x 251 x 80 mm              |
| B2-4825-100 | 483 x 251 x 110 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-4825  | 485 x 256 x 5 mm               |

Quarter Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-4025-30  | 403 x 251 x 40 mm              |
| B2-4025-50  | 403 x 251 x 60 mm              |
| B2-4025-70  | 403 x 251 x 80 mm              |
| B2-4025-100 | 403 x 251 x 110 mm             |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-4025  | 405 x 256 x 5 mm               |

Dental Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2717-40  | 273 x 172 x 50 mm              |
| B2-2717-100 | 273 x 172 x 110 mm             |
| -           | -                              |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-2717  | 275 x 174 x 5 mm               |

Mini Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-2712-40  | 260 x 120 x 50 mm              |
| B2-2712-100 | 260 x 120 x 110 mm             |
| -           | -                              |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-2712  | 262 x 122 x 5 mm               |

Dental 1/2



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B2-1714-40  | 165 x 130 x 45 mm              |
| B2-1714-100 | 165 x 130 x 110 mm             |
| -           | -                              |
| -           | -                              |

| Lid Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B2-1714  | 167 x 130 x 5 mm               |



Half Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-2525-50  | 245 x 255 x 50 mm              |
| B3-2525-70  | 245 x 255 x 70 mm              |
| B3-2525-100 | 245 x 255 x 100 mm             |
| Lid Code    | Outside Dimension<br>L x W x H |
| B3-2525     | 215 x 240 x 5 mm               |

Long Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-6625-50  | 665 x 255 x 50 mm              |
| B3-6625-70  | 665 x 255 x 70 mm              |
| B3-6625-100 | 665 x 255 x 100 mm             |
| Lid Code    | Outside Dimension<br>L x W x H |
| B3-6625     | 635 x 240 x 5 mm               |

Wide Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-5333-50  | 533 x 337 x 50 mm              |
| B3-5333-70  | 533 x 337 x 70 mm              |
| B3-5333-100 | 533 x 337 x 100 mm             |
| Lid Code    | Outside Dimension<br>L x W x H |
| B3-5333     | 503 x 322 x 5 mm               |

Fine Mesh Baskets

Full Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-5425-50  | 540 x 255 x 50 mm              |
| B3-5425-70  | 540 x 255 x 70 mm              |
| B3-5425-100 | 540 x 255 x 100 mm             |
| Lid Code    | Outside Dimension<br>L x W x H |
| B3-5425     | 510 x 240 x 5 mm               |

Din Size





| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-4825-50  | 485 x 255 x 50 mm              |
| B3-4825-70  | 485 x 255 x 70 mm              |
| B3-4825-100 | 485 x 255 x 100 mm             |
| Lid Code    | Outside Dimension<br>L x W x H |
| B3-4825     | 455 x 240 x 5 mm               |

Quarter Size



| Box Code    | Outside Dimension<br>L x W x H |
|-------------|--------------------------------|
| B3-4025-50  | 405 x 255 x 50 mm              |
| B3-4025-70  | 405 x 255 x 70 mm              |
| B3-4025-100 | 405 x 255 x 100 mm             |
| Lid Code    | Outside Dimension<br>L x W x H |
| B3-4025     | 375 x 240 x 5 mm               |




| Box Code | Outside Dimension<br>L x W x H |
|----------|--------------------------------|
| B4-33-30 | 30 x 30 x 30 mm                |
| B4-44-20 | 40 x 40 x 20 mm                |
| B4-66-20 | 60 x 60 x 20 mm                |
| B4-66-40 | 60 x 60 x 40 mm                |
| B4-84-40 | 80 x 40 x 40 mm                |

| Box Code   | Outside Dimension<br>L x W x H |
|------------|--------------------------------|
| B4-88-40   | 80 x 80 x 40 mm                |
| B4-88-60   | 80 x 80 x 60 mm                |
| B4-1212-30 | 120 x 120 x 30 mm              |
| B4-1212-60 | 120 x 120 x 60 mm              |
| B4-2424-60 | 240 x 240 x 60 mm              |

Scope Baskets



Perforated Sheet  
Small Size Baskets



 Bottom With Lid  
W/O Feet  
Wire:  $\varnothing$  1 mm  
Mesh: 4x4 mm

| Code       | Dimension<br>L x W x H | Code      | Dimension<br>L x W x H |
|------------|------------------------|-----------|------------------------|
| B2-2718-32 | 275 x 180 x 32 mm      | B2-188-32 | 188 x 86 x 32 mm       |
| B2-2718-42 | 275 x 180 x 42 mm      | B2-144-28 | 146 x 46 x 28 mm       |
| B2-247-30  | 240 x 76 x 30 mm       | B2-84-28  | 86 x 46 x 28 mm        |
| B2-1913-34 | 191 x 130 x 34 mm      |           |                        |

STU - SPRI Baskets

| Wire Basket<br>Single Scope                                                                                                                    | Wire Basket<br>Double Scope                                                                                                                    | Perforated Basket<br>Single Scope                                                                                                              | Perforated Basket<br>Double Scope                                                                                                                  |
|------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Code<br>B5-298-55-S<br><br>Dimension<br>L x W x H<br>290 x 80 x 52 mm<br><br>Code<br>B5-468-55-S<br>Dimension<br>L x W x H<br>460 x 80 x 52 mm | Code<br>B5-298-55-D<br><br>Dimension<br>L x W x H<br>290 x 80 x 52 mm<br><br>Code<br>B5-468-55-D<br>Dimension<br>L x W x H<br>460 x 80 x 52 mm | Code<br>B5-278-44-S<br><br>Dimension<br>L x W x H<br>272 x 83 x 44 mm<br><br>Code<br>B5-458-44-S<br>Dimension<br>L x W x H<br>452 x 83 x 44 mm | Code<br>B5-2710-44-D<br><br>Dimension<br>L x W x H<br>272 x 103 x 44 mm<br><br>Code<br>B5-4510-44-D<br>Dimension<br>L x W x H<br>452 x 103 x 44 mm |

| STU Baskets | Box Code    | Outside Dimension L x W x H | Details                                                             |
|-------------|-------------|-----------------------------|---------------------------------------------------------------------|
|             | B7-5728-135 | 575 x 280 x 135 mm          | Upper Frame: $\varnothing$ 6 mm                                     |
|             | B7-5728-265 | 575 x 280 x 265 mm          | Body Wires: $\varnothing$ 2 & 3 mm<br>Base Frame $\varnothing$ 5 mm |

| SPRI Baskets | Box Code    | Outside Dimension L x W x H | Details                        |
|--------------|-------------|-----------------------------|--------------------------------|
|              | B8-6030-40  | 600x300x40 mm               | Main Frame: $\varnothing$ 6 mm |
|              | B8-6030-100 | 600x300x100 mm              | Body Wire: $\varnothing$ 4 mm  |
|              | B8-6030-150 | 600x300x150 mm              | Base Frame $\varnothing$ 5 mm  |
|              | B8-6040-100 | 600x400x100 mm              | Single Frame Stackable         |
|              | B8-6040-150 | 600x400x150 mm              |                                |
|              | B8-6040-200 | 600x400x200 mm              |                                |



| Paper Filters                                                                                     |                                                             |
|---------------------------------------------------------------------------------------------------|-------------------------------------------------------------|
|                                                                                                   |                                                             |
| Code A1-10-1<br>For: 1/1, 3/4, 1/2, Long, Wide<br>Size: ø 185 mm                                  | Code A1-11-1<br>For: Small Size Containers<br>Size: ø 88 mm |
| Material 73g/m² Medical Grade Paper<br>Package 100 pcs<br>Sterilization Steam<br>Cycle Single Use |                                                             |

| PTFE Teflon Filters                                                                             |                                                                                            |
|-------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
|                                                                                                 |                                                                                            |
| Code A1-20-1<br>For: 1/1, 3/4, 1/2, Long, Wide<br>Size: ø 190 mm<br>Cycle: 2000 Cycles Reusable | Code A1-21-1<br>For: Small Size Containers<br>Size: ø 88 mm<br>Cycle: 2000 Cycles Reusable |
| Code A1-20-2<br>For: 1/1, 3/4, 1/2, Long, Wide<br>Size: ø 190 mm<br>Cycle: 1200 Cycles Reusable | Code A1-21-2<br>For: Small Size Containers<br>Size: ø 88 mm<br>Cycle: 1200 Cycles Reusable |

| Polymer Security Seals                                                                   |                                                                                         |
|------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|                                                                                          |                                                                                         |
| Code A2-10-1<br>Package 100 pcs<br>Sterilization Steam<br>Cycle Single Use<br>Color Blue | Code A2-20-1<br>Package 100 pcs<br>Sterilization Steam<br>Cycle Single Use<br>Color Red |

| Polymer Security Seals With Indicator                                                    |                                                                                         |
|------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|                                                                                          |                                                                                         |
| Code A2-10-2<br>Package 100 pcs<br>Sterilization Steam<br>Cycle Single Use<br>Color Blue | Code A2-20-2<br>Package 100 pcs<br>Sterilization Steam<br>Cycle Single Use<br>Color Red |



| Identification Labels                                                                               |         |       |         |
|-----------------------------------------------------------------------------------------------------|---------|-------|---------|
| Color                                                                                               | Code    | Color | Code    |
|                                                                                                     | A5-10-1 |       | A5-10-5 |
|                                                                                                     | A5-10-2 |       | A5-10-6 |
|                                                                                                     | A5-10-3 |       | A5-10-7 |
|                                                                                                     | A5-10-4 |       | A5-10-8 |
| Material 1,5 mm Aluminum<br>Package 16 x 25 mm<br>Sterilization Steam<br>Marking Maximum 24 Letters |         |       |         |

| Paper Labels                                                                 |                                |
|------------------------------------------------------------------------------|--------------------------------|
|                                                                              |                                |
| Code A3-10-1<br>with indicator                                               | Code A3-10-2<br>with indicator |
| Material Paper<br>Size Big & Small<br>Package 100 pcs<br>Sterilization Steam |                                |

| Filter Retainers                               |                                            |
|------------------------------------------------|--------------------------------------------|
|                                                |                                            |
| Code A6-10-1<br>For: 1/1, 3/4, 1/2, Long, Wide | Code A6-10-2<br>For: Small Size Containers |
| Material Stainless Steel<br>Package 1 piece    |                                            |

| Silicone Mats & Nets |                                   |          |                                                                                               |
|----------------------|-----------------------------------|----------|-----------------------------------------------------------------------------------------------|
|                      | Code A7-5425<br>Size 500 X 235 mm |          | Code A4-5425<br>Size 540 x 252 mm                                                             |
|                      | A7-4025 370 X 235 mm              |          | A4-4025 420 x 252 mm                                                                          |
|                      | A7-2525 225 x 235 mm              |          | A4-2525 254 x 253 mm                                                                          |
|                      | A7-2712 270 x 125 mm              |          |                                                                                               |
|                      | h: 14 mm<br>h: 5 mm<br>h: 2 mm    | h: 16 mm | Material Medical Grade Silicone<br>Package 1 Piece<br>Sterilization Steam<br>Color Light Blue |

| Silicone Gasket                                                                                                                |                                |                                      |                                |                            |                            |
|--------------------------------------------------------------------------------------------------------------------------------|--------------------------------|--------------------------------------|--------------------------------|----------------------------|----------------------------|
|                                                                                                                                | Code A11-10-1<br>Size 1/1 Full | Code A11-20-13/4<br>Size 1/4 Quarter | Code A11-30-1<br>Size 1/2 Half | Code A11-40-1<br>Size Long | Code A11-50-1<br>Size Wide |
|                                                                                                                                | A11-60-1 Dental                | A11-70-1 Mini                        | A11-80-11/2 Dental             | A11-90-1 Flat              |                            |
| Material Medical Grade Sealing Silicone<br>Package 1 Piece<br>Shape D shape, Jointless<br>Feature Place & Go, No Glue Required |                                |                                      |                                |                            |                            |



Code T1-10-1  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece  
Dimensions: 600 x 860 x 900 mm (W x L x H)



Code T1-10-2  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece  
Dimensions: 550 x 860 x 900 mm (W x L x H)



Code T1-10-5  
Material : Stainless Steel  
Feet : 4 Adjustable Leveling Feet  
Package : 1 Piece  
Dimensions: 600 x 1400 x 850 mm (W x L x H)




Code T1-10-6  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece  
Dimensions: 700 x 360 x 1350 mm (W x L x H)



Code T1-10-3  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece  
Dimensions: 600 x 850 x 900 mm (W x L x H)



Code T1-10-4  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece  
Dimensions: 550 x 860 x 900 mm (W x L x H)



Code T1-10-7  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece  
Dimensions: 700 x 1080 x 1350 mm (W x L x H)



Code T1-10-8  
Material : Stainless Steel  
Feet : 4 wheels - 2 braked  
Package : 1 Piece - including 3 baskets  
Dimensions: 400 x 600 x 1700 (W x L x H)